



The Laboratory People

# camlab

Camlab House, Norman Way Industrial Estate  
Over, Cambridge, CB24 5WE, United Kingdom

Tel: +44 (0)1954 233 100 Fax: +44 (0)1954 233 101  
Email: sales@camlab.co.uk Web: www.camlab.co.uk



# LABCOLD™

medical & scientific refrigeration



## RAVF1825 Advanced Laboratory Freezer

| Part Number                 | RAVF1825           |
|-----------------------------|--------------------|
| Capacity (Litres)           | 505                |
| Style                       | Freestanding       |
| Exterior Dimensions (HxWxD) | 1702 x 770 x 725mm |
| Interior Dimensions (HxWxD) | 1490 x 645 x 545mm |
| Weight                      | 86kg               |
| Operating Temperature       | -20°C              |
| Interior Construction       | White Enamel Steel |
| Exterior Construction       | White Enamel Steel |
| Fixed Shelves               | 6                  |
| Door Type                   | Solid              |
| Door Hinging                | Right              |
| Reverse Hinged              | Optional           |
| Defrost                     | Manual             |
| Access Port                 | Optional           |
| Open Door Alarm             | X                  |
| High Temperature Alarm      | ✓                  |
| Low Temperature Alarm       | ✓                  |
| Remote Alarm Contacts       | Optional           |
| Digital Temperature Display | ✓                  |
| Door Lock                   | ✓                  |
| Feet                        | 4                  |
| Energy Consumption          | 4.1 kWh/24h        |
| Eco Friendly Refrigerant    | ✓                  |

Labcold Advanced freezers are constructed from powder coated mild steel with an enamel steel interior so are built to last. Ideal for any application which requires the storage of a large quantity of product within tight temperature tolerances. This freezer is even manual defrost, putting you completely in control of the temperature of your samples.

**Looking for  
temperature  
mapping or  
calibration?  
Call the experts on  
01256 705 570.**

Made to order at our UK factory, this freezer is specifically designed to withstand the rigours of professional use, and is backed by a comprehensive one year parts and labour warranty.

Fitted with a digital temperature display, you can be sure with a quick glance that the freezer is operating at the correct temperature. For extra peace of mind this freezer also has an audio/visual minimum and maximum temperature alarm and can have remote alarm contacts fitted at time of manufacture. For extra security this freezer comes complete with a door lock and two keys.

01256 705 570

www.labcold.com

