

# SK-O180-Pro/SK-L180-Pro Orbital and Linear Shakers




## Features

- Compact size
- LCD display
- Electronic Timer 1-1199 mins
- Free ShakerPC software

The Camlab Choice SK-0180/SK-L180-Pro Orbital and Linear shakers provide smooth and quiet orbital or linear mixing for a wide range of solutions, culture tubes and flasks. They are compact in size and have an electronic timer for consistent shaking between batches. The shaker also has a continuous run mode. Free software ShakerPC is available to control the unit through an RS232 interface.

## Specifications

Part Number	Linear Shaker	SK-L180-Pro	1189970
	Orbital Shaker	SK-O180-Pro	1189972
Speed range	Linear Shaker	SK-L180-Pro	100-350 rpm
	Orbital Shaker	SK-O180-Pro	100-800 rpm
Orbital diameter			4mm
Frequency			50/60Hz
Max load			2.5Kg
Dimensions			340 x 300 x 100
Weight			8.1Kg
Protection class			IP21
Power			30W



## Accessories

### Cat No

**1191820 Universal attachment with 3 bars for SK-O180/SK-L180Pro Orbital/Linear Shaker**

