



Model:  
**ST-155**

Typology:  
**STAND FOR MODULAR STEREO MICROSCOPES**

**Description:**  
*Large plain stand with head holder and focusing mechanism. Provided with incident and transmitted LED illuminator.*

<b>Illumination</b>	<p>White LED incident and transmitted illumination system, with independent intensity controls using knobs at both sides of the base.          LED power 3W, comparable to an halogen bulb 30W.          Color temperature: 6300K          LED average life time approx. 50000h.          Voltage: 110/240Vac, 50/60Hz, 0,45 A; Fusibile: T2A 250V          Max. power required: 7W.</p>
<b>Focusing</b>	<p>By means of a rack and pinion controlled by a pair of knobs placed on both sides of the head holder.          The tension of focusing can be adjusted.          The head holder can slide along the fixed arm.</p>
<b>Dimensions</b>	<p>Base: 330 x 285 mm          Fixed arm: height 305 mm.</p>