




Meeting Your Needs in Every Weigh



Sartorius Talent Series

**Portable,
Dependable,
Accurate**

Sartorius Talent Series

A “Star Performer” in Basic Weighing

Rugged and compact, the Talent series from Sartorius answers the need for a portable, lightweight, basic balance that can go anywhere and give accurate results everytime.

Nineteen models to choose from is just the start

The Talent series offers 3 analytical models and 16 toploaders. Weighing capacities go from 60g to 12kg. Accuracy/readability goes from 0.0001g to 1g.

Easy to use, portable with dependable results

“Plug and play” and the Talent is ready to go - in the lab or in the field. Rugged construction with a stainless steel pan, and sealed overlay assures years of usage.

Model	Weighing capacity (g)	Readability (g)	Pan size (mm/in)	Response time average (s)	Repeatability (< ± g)	Linearity (< ± g)	Fax Doc. No.	Battery Operable (8AA Batteries)
Analyticals								
TE64	60	0.0001	∅ 90 / 3.5	3	0.0001	0.0002	10072	No
TE124S	120	0.0001	∅ 90 / 3.5	3	0.0001	0.0002	10071	No
TE214S	210	0.0001	∅ 90 / 3.5	3	0.0001	0.0002	10070	No
Toploaders								
TE153S	150	0.001	∅ 100 / 3.9	2.5	0.0015	0.003	10075	No
TE153S-DS *	150	0.001	∅ 100 / 3.9	2.5	0.0015	0.003	10076	No
TE313S	310	0.001	∅ 100 / 3.9	2.5	0.001	0.002	10073	No
TE313S-DS *	310	0.001	∅ 100 / 3.9	2.5	0.001	0.002	10074	No
TE212	210	0.01	∅ 116 / 4.5	2	0.01	0.02	10081	Yes
TE412	410	0.01	∅ 116 / 4.5	2	0.01	0.02	10080	Yes
TE612	610	0.01	∅ 116 / 4.5	2	0.01	0.02	10079	Yes
TE1502S	1,500	0.01	174 x 143 / 6.9 x 5.6	2.5	0.015	0.03	10078	No
TE3102S	3,100	0.01	174 x 143 / 6.9 x 5.6	2.5	0.01	0.02	10077	No
TE601	610	0.1	174 x 143 / 6.9 x 5.6	1.5	0.1	0.2	10085	Yes
TE2101	2,100	0.1	174 x 143 / 6.9 x 5.6	1.5	0.1	0.2	10084	Yes
TE4101	4,100	0.1	174 x 143 / 6.9 x 5.6	2	0.1	0.2	10083	Yes
TE6101	6,100	0.1	174 x 143 / 6.9 x 5.6	2	0.1	0.2	10082	Yes
TE4100	4,100	1	174 x 143 / 6.9 x 5.6	1.5	1	2	10088	Yes
TE6100	6,100	1	174 x 143 / 6.9 x 5.6	1.5	1	2	10087	Yes
TE12000	12,000	1	174 x 143 / 6.9 x 5.6	1.5	1	2	10086	Yes

* Includes standard analytical draft chamber. Call 800-551-3774 and enter the fax number for more information.

Advanced Features

- Monolith DLS weighing system for fast stable results
- Standard RS-232 adapter connects to data printer, PC or to a remote display
- Rugged housings and stainless steel weighing pans

Standard Applications

- Weighing in percent
- Net-total formulation
- Animal weighing (averaging)
- Counting
- Toggling between two weight units

Selected Accessories

- YDP03-OCE
Statistical data printer
- YSC01L
SartoConnect data transfer software
- YRB08Z
Rechargeable battery pack

Sartorius Corp.
131 Heartland Blvd.
Edgewood, NY 11717
800-635-2906
www.sartorius.com

Sartorius Inc.
2179 Dunwin Drive #4&5
Mississauga, Ontario
Canada L5L 1X2
800-668-4234



sartorius

Meeting Your Needs in Every Weigh