



# Flasks



## Flasks, conical, wide neck

**PYREX®**

- Complies with ISO 24450
- Erlenmeyer pattern
- Manufactured from Pyrex® borosilicate glass
- Excellent thermal performance
- High resistance to chemical attack
- White enamel graduations and marking spot

Product Code	Capacity (ml)	Body O.D. mm	Height mm	Neck I.D. mm	Pack qty.
1140/02M	50	48	85	30	10
1140/04D	100	64	110	30	10
1140/08D	250	85	140	46	10
1140/10D	500	105	175	46	10
1140/14D	1000	131	220	46	10
1140/18D	2000	153	275	66	10

**i** 2000ml size is supplementary to ISO range.



## Flasks, conical, wide neck, heavy duty

**PYREX®**

- Mechanically stronger than conventional Erlenmeyers
- Reinforced rims for extra strength and resistance to impact fracture
- Manufactured from high quality Pyrex® borosilicate glass
- White graduations and marking spot

Product Code	Capacity (ml)	Height mm	Body O.D. mm	Neck I.D. mm	B.S. stopper size	Pack qty.
5100-125	125	109	66	27	27	12
5100-250	250	134	81	33	33	12
5100-500	500	170	97	43	43	6
5100-1L	1000	219	122	48	48	6
5100-2L	2000	279	156	58	58	1

