



Funnels, dropping, Rotaflor[®] stopcock, stem with cone



- Complies to ISO 4800
- Funnel stem with ground-glass cone for easy incorporation into the sockets of Quickfit[®] flasks etc.
- Fitted with GP Rotaflor[®] stopcock with PTFE key
 - No lubrication required
 - No contaminating greases
 - No seizure with caustic solutions

Product Code	Nom. capacity ml	Stopcock bore mm	Cone size	Socket size	Pack qty.
D1/11GP	50	3.0	14/23	14/23	1
D1/21GP	100	3.0	14/23	14/23	1
D1/22GP	100	3.0	19/26	19/26	1
D1/32GP	250	3.0	19/26	24/29	1
D1/42GP	500	6.0	19/26	24/29	1
D1/62GP	1000	6.0	19/26	29/32	1



Funnels, dropping, Rotaflor[®] stopcock, open top, with cone



- Complies with ISO 4800
- With stem inside the cone
- Enables easy fitting to the top of Quickfit[®] chromatography columns etc where a protruding stem would be undesirable
- GP Rotaflor[®] stopcock with PTFE key for grease-free performance
- With open top instead of socket for easy addition of solutions
- Economical price

Product Code	Nom. capacity ml	Top I.D. mm	Stopcock bore mm	Cone size	Pack qty.
D3/12	50	28	3.0	19/26	1
D3/22	100	39	3.0	19/26	1
D3/32	250	47	3.0	19/26	1

