



## Funnels, dropping, graduated, glass key

**PYREX®**

- Complies with ISO 4800 type 4
- Manufactured from Pyrex® borosilicate glass
- With double graduated scale showing approximate volume contained/dispensed
- Graduations in white enamel
- With interchangeable glass key
- Supplied with glass stopper

Product Code	Capacity (ml)	Grads. ml	Stopcock bore mm	Stopper size	Pack qty.
3590/01	50	1	1.5	19/26	1
3590/02	100	2	2.5	19/26	1
3590/04	250	5	2.5	24/29	1
3590/06	500	10	4.0	24/29	1
3590/08	1000	20	4.0	29/32	1



Spare glass stoppers are detailed on page 247.



## Funnels, dropping, graduated, glass key, with cone

**PYREX®**

- Complies with ISO 4800
- With interchangeable glass key
- Double graduated scale showing volume of liquid contained/delivered
- Stem inside cone for easy fitting to Quickfit® assemblies
- Graduations in white enamel

Product Code	Capacity (ml)	Grads. ml	Stopcock bore mm	Cone and socket size	Pack qty.
443334	50	1	1.5	19/26	1
443341	100	2	1.5	19/26	1
443348	250	5	2.5	29/32	1
443355	500	10	4.0	29/32	1
443361	1000	20	4.0	29/32	1