



Funnels, dropping, glass key, stem with cone



- Complies with ISO 4800
- Funnel stem with ground-glass cone for easy incorporation into Quickfit® flasks and assemblies
- Stopcock has interchangeable glass key

Product Code	Nom. capacity ml	Stopcock bore mm	Cone size	Socket size	Pack qty.
D1/00	20	1.5	10/19	10/19	1
D1/11	50	1.5	14/23	14/23	1
D1/21	100	2.5	14/23	14/23	1
D1/22	100	2.5	19/26	19/26	1
D1/32	250	2.5	19/26	24/29	1
D1/42	500	4.0	19/26	24/29	1
D1/62	1000	4.0	19/26	29/32	1



i Stoppers for use with Quickfit® dropping funnels are detailed on page 247.

Funnels, pressure equalising, pear shape, Rotaflo® stopcock



- Complies with ISO 4800
- For use when liquids have to be introduced into vessels under vacuum
- Can also be used with Soxhlets as solvent traps to facilitate recovery at the end of an extraction
- With GP Rotaflo® stopcocks having chemically resistant PTFE keys

Product Code	Nom. capacity ml	Stopcock bore mm	Cone size	Socket size	Pack qty.
DE100/22	100	3.0	19/26	19/26	1
DE100/32	100	3.0	19/26	24/29	1
DE250/32	250	3.0	19/26	24/29	1
DE250/33	250	3.0	24/29	24/29	1
DE250/44	250	3.0	29/32	29/32	1
DE250/55	250	3.0	34/35	34/35	1

