



Funnels, pressure equalising, cylindrical, graduated, PTFE key **PYREX®**

- Complies with ISO 4800
- With interchangeable PTFE key for grease free performance
- Double graduated scale showing volume if liquid contained/delivered
- For use when liquids have to be added to systems under vacuum
- Stem inside cone for easy fitting to Quickfit® assemblies
- Graduations in white enamel
- Supplied with hollow glass stopper

Product Code	Capacity (ml)	Grads. ml	Stopcock bore mm	Cone and socket size	Pack qty.
443533	50	1	1.5	19/26	1
443540	100	2	1.5	19/26	1
443547	250	5	2.5	29/32	1
443554	500	10	4.0	29/32	1
443560	1000	20	4.0	29/32	1



Funnel, pressure equalising, cylindrical, graduated, glass key **PYREX®**

- Complies with ISO 4800
- With interchangeable glass key
- Double graduated scale showing volume if liquid contained/delivered
- For use when liquids have to be added to systems under vacuum
- Stem inside cone for easy fitting to Quickfit® assemblies
- Graduations in white enamel
- Supplied with hollow glass stopper

Product Code	Capacity (ml)	Grads. ml	Stopcock bore mm	Cone and socket size	Pack qty.
443534	50	1	1.5	19/26	1
443541	100	2	1.5	19/26	1
443548	250	5	2.5	29/32	1
443555	500	10	4.0	29/32	1
443561	1000	20	4.0	29/32	1