



Funnels, Industrial, HDPE



- Rigid industrial funnels in polyethylene
- Unbreakable in normal use
- Feature straight sided upper position to prevent spillage

Item Code	Top I.D. mm	Stem O.D. at tip mm	Overall height mm	Pack qty.
FWC102	80	9	80	10
FWC104	100	11	102	10
FWC106	120	12	121	10
FWC108	140	13	141	5
FWC110	180	17	181	5
FWC112	230	21	226	5
FWC114	245	24	239	1
FWC116	290	27	277	1
FWC118	360	36	351	1



Funnels, powder, polypropylene



- Lightweight polypropylene powder funnel
- Unbreakable in normal use
- Parallel stem minimises 'bridging' of powder
- The two largest sizes have vertical sided tops

Item Code	Top I.D. mm	Stem O.D. mm	Overall height mm	Pack qty.
FWC202	65	18	69	5
FWC204	80	22	74	5
FWC206	100	26	90	5
FWC208	120	30	106	5
FWC210	150	36	135	5
FWC212	290	35	208	1
FWC214	420	48	270	1