



## Beakers, straight sided, printed graduations, polypropylene



- Popular classic Azlon beakers
- Ideal for general laboratory use - straight stable shape
- Polypropylene material - good chemical resistance
- Can withstand temperatures up to 121°C and therefore autoclavable
- Graduations are clearly printed in indelible solvent resistant ink

Item Code	Capacity ml	Grads. ml	O.D. mm	Height mm	Pack qty.
BDA208P	10	1.0	32	32	10
BDA212P	25	2.0	42	40	10
BDA218P	50	5.0	51	51	10
BDA224P	100	10.0	62	62	10
BDA230P	250	25.0	82	82	10
BDA234P	400	50.0	94	95	10
BDA236P	500	50.0	100	101	10
BDA238P	600	50.0	106	107	5
BDA240P	1000	100.0	127	126	5
BDA242P	2000	200.0	158	160	1



## Beakers, straight sided, printed graduations, polypropylene, sets

Item Code	Description	Pack qty.
BDA218MES	Set of 3 beakers (50, 100 & 250ml)	1
BDA1198	Set of 3 beakers (500, 1000 & 2000ml)	1
BDA1199	Set Of 5 beakers (10, 25, 50, 100 & 250ml)	1
BDA1200	Set of 8 beakers (10, 25, 50, 100, 250, 500, 1000 & 2000ml)	1
BDA300P	Set of 10 beakers (One of each size)	1