



## Beakers, straight sided, printed graduations, polypropylene, with dust covers



- Popular classic Azlon beaker supplied with tight fitting dust cover lid
- Protects against airborne contamination and accidental splashing or spillage

Item Code	Capacity ml	Grads. ml	O.D. mm	Height mm	Pack qty.
BDA730P	250	25.0	82	82	5
BDA738P	600	50.0	106	107	5
BDA740P	1000	100.0	127	126	1
BDA742P	2000	200.0	158	160	1



*BDA740P and BDA742P have non-autoclavable polystyrene lids.*



## Beakers, tapered, printed graduations, polypropylene



- Tapered stacking style - space saving and easy to handle
- Comply with ISO 7056
- High clarity polypropylene - autoclavable at 121°C with good chemical resistance
- Crisp, blue, easy to read indelible graduations
- Precise non drip pouring spouts

Item Code	Capacity ml	Grads. ml	Rim diam. mm	Height mm	Pack qty.
BPB0025P	25	5.0	41.5	51	10
BPB0050P	50	10.0	50.5	61	10
BPB0100P	100	20.0	61.5	73	10
BPB0250P	250	50.0	80.5	96	10
BPB0400P	400	50.0	92.5	111	10
BPB0500P	500	50.0	102.5	121	10
BPB0600P	600	50.0	104.0	126	5
BPB1000P	1000	100.0	120.0	148	5
BPB2000P	2000	200.0	163.0	188	1
BPB3000P	3000	200.0	172.0	214	1
BPB5000P	5000	500.0	230.0	250	1