

RANGE 11.A

VENTILATED SAFETY CABINETS FOR ACIDS AND BASES



CONFORMITY

- Large warning labels in compliance with ISO 3864, 92/58/CEE European directive and NF X08.003 norm.
- Electric fan in compliance with CE regulations.

ADVANTAGES

- Strictly separated acid and base storage.

PASSIVE SECURITY

- Built in steel 15/10th.
- Double wall doors in steel 15/10th.
- White epoxy paint RAL 9010.
- «Acid» and «Base» large warning labels.
- Retention shelves in steel adjustable in height.




- Additional retention shelves in PVC for «Acid» compartments.
- Removable base through in steel at the bottom of each compartment (except «Acid» compartments: base through in PVC).

ACTIVE SECURITY

- High and low vents Ø 100 mm (for possible external connection).
- Handling doors with key-locking system.
- Electric fan: 220/240 volts - 50 Hz - 40 W - 0,19 A - 54 dB(A) - Flow: 156 m³/h - Protection: IP44.

OPTIONS

- Industrial absorbing mat (ref. JF10D).
- Flexible tube (ref. KL100).
- Clamp collar (ref. CDS100).

Ref.	Description	Dimensions H x W x D (mm)	Compartment capacity		Nb of shelves	Nb of base through	Fan	Compartment type			W. (kg)
			Stockage (L)	Rétention (L)							
AB75	Tall cabinet with fan	Ext. 2080x560x560 Int. (2x) 890x530x470	110 (2x 55)	70 (2x 35)	4	2	156 m³/h	x	x	-	78
AB75SC	Tall cabinet without fan	Ext. 1830x560x560 Int. (2x) 890x530x470					-	x	x	-	68
AS75SC	Tall cabinet without fan						-	x	-	x	68