

# premium pipette tips

## Qualitips®

From the 10 µL ultra-microtip to the 10 mL macrotip, the Qualitips® range offers a wide selection of autoclavable superior quality pipette tips covering most research and routine applications. More than just compatible with the Socorex micropipettes, the lot by lot quality controlled Qualitips® fit also a variety of other brands.

### Qualitips® without filter, up to 1000 µL

- Made from virgin, metal free polypropylene
- Precision molded for air tight fit
- Autoclavable
- Supplied in bags and autoclavable racks with lid

Cat. No.	Vol.	Description	Packaging
----------	------	-------------	-----------


Fit Socorex®, Gilson®



Length : 31.0 mm

309.0010B	10 µL	Natural ultra-microtip	1x 1000/bag
309.0010R	10 µL	Natural ultra-microtip	10x 96/rack


Fit Socorex®, Eppendorf®



Length : 42.3 mm

302.0020B	10 µL	Natural microtip	1x 1000/bag
302.0020R	10 µL	Natural microtip	10x 96/rack


Fit Socorex®, Eppendorf®



Length : 70.9 mm

302.0010GR	10 µL	Natural gel loading tip	4x 204/rack
------------	-------	-------------------------	-------------

Fit Socorex®, Eppendorf®, Gilson®



Length : 49.4 mm

307.0200B	200 µL	Natural universal tip	2x 500/bag
307.0200R	200 µL	Natural universal tip	10x 96/rack
327.0200B	200 µL	Yellow universal tip	2x 500/bag
327.0200R	200 µL	Yellow universal tip	10x 96/rack

Fit Socorex®, Eppendorf®, Gilson®



Length : 51.1 mm

308.0200B	200 µL	Natural tip	1x 1000/bag
308.0200R	200 µL	Natural tip	10x 96/rack

Fit Socorex®, Gilson®



Length : 49.7 mm

309.0200B	200 µL	Natural tip	2x 500/bag
309.0200R	200 µL	Natural tip	10x 96/rack
329.0200B	200 µL	Yellow tip	2x 500/bag
329.0200R	200 µL	Yellow tip	10x 96/rack

Cat. No.	Vol.	Description	Packaging
----------	------	-------------	-----------

Fit Socorex®, Eppendorf®



Length : 49.0 mm

328.0200B	200 µL	Yellow tip	2x 500/bag
328.0200R	200 µL	Yellow tip	10x 96/rack


Fit Socorex®, Eppendorf®, Gilson®



Length : 69.0 mm

308.0200GR	200 µL	Natural gel loading tip	4x 204/rack
------------	--------	-------------------------	-------------


Fit Socorex®, Eppendorf®, Gilson®



Length : 58.8 mm

308.0350B	350 µL	Natural tip	1x 1000/bag
308.0350R	350 µL	Natural tip	10x 96/rack

Fit Socorex®, Eppendorf®, Gilson®



Length : 79.1 mm

307.1000B	1000 µL	Natural universal tip	4x 250/bag
307.1000R	1000 µL	Natural universal tip	10x 60/rack


Fit Socorex®, Gilson®



Length : 72.0 mm

309.1000B	1000 µL	Natural tip	4x 250/bag
309.1000R	1000 µL	Natural tip	8x 60/rack
319.1000B	1000 µL	Blue tip	4x 250/bag
319.1000R	1000 µL	Blue tip	8x 60/rack

Fit Socorex®, Eppendorf®






Length : 69.5 mm


318.1000B	1000 µL	Blue tip	4x 250/bag
318.1000R	1000 µL	Blue tip	8x 60/rack

## Qualitips® without filter, up to 10 mL

- Natural macro tips made from virgin polypropylene
- Dye and metal free
- Autoclavable
- Supplied in bags and autoclavable racks with lid

	Cat. No.	Vol.	Description	Packaging
<p><i>Fit Socorex®, Eppendorf®</i></p>  <p>Length : 115.5 mm</p>	312.02	2 mL	Natural macrotip	2x 250/bag
<p><i>Fit Socorex®, Eppendorf®</i></p>  <p>Length : 123.2 mm</p>	312.05B 312.05R	5 mL 5 mL	Natural macrotip Natural macrotip	2x 250/bag 2x 50/rack
<p><i>Fit Socorex®, FinnpiPETTE®</i></p>  <p>Length : 150.0 mm</p>	312.10	10 mL	Natural macrotip	3x 100/bag


## Glass Pasteur pipette

 <p>Length : 150.0 mm</p>	313.02	2 mL	Pasteur pipette, glass	250/pack
---	--------	------	------------------------	----------

## Nozzle protection filter for macropipettes

- Cellulose filter for added protection against liquid fill and contamination

Cat. No.	Volume	Description	Packaging
322.02	2 mL	For pipettes Calibra® 832 and Acura® analog 831	1 x 250/bag
322.05	5 mL	For pipettes Acura® 835 and 935	1 x 250/bag
322.10	10 mL	For pipettes Calibra® 832, Acura® 835 and 935	1 x 100/bag



## Autoclaving boxes and all purpose reagent reservoirs


- Well adapted for single and multichannel pipettes
- Various shapes in different configurations





Cat. No.	Description	Packaging
<b>Autoclaving boxes with lid (supplied empty)</b>		
302.0020ER	For Qualitips 302.0020, 10 µL	10/pack
307.0200ER	With hinged lid for Qualitips 307.00200, 200 µL	10/pack
328.0200ER	For Qualitips 328.0200, 200 µL	10/pack
307.1000ER	With hinged lid for Qualitips 307.1000, 1000 µL	10/pack
319.1000ER	For Qualitips 319.1000, 1000 µL	8/pack
312.10ER	For Qualitips 312.10, 10 mL	4/pack
<b>Reagent reservoirs</b>		
330.01	V-shape, PP, 70 mL, autoclavable	20/pack
330.01L	V-shape with hinged lid, PP, 175 mL, autoclavable	20/pack
330.08.9	8-channel, PS, 40 mL, sterile, single wrapped	10 x 1/pack
330.12.9	12-channel, PS, 40 mL, sterile, single wrapped	10 x 1/pack


## Qualitips® with filter up to 1000 µL


- Made from premium dye and metal free polypropylene
- Certified DNase + RNase free. sterile tips pyrogen free
- Specially designed to prevent aerosol contamination
- High density PE filter allows complete, uncontaminated sample retrieval
- Lot by lot quality control
- Supplied in racks with lid or bags, sterile or autoclavable


Cat. No.	Vol.	Description	Packaging
<i>Fit Socorex®, Gilson®</i>			
			
Length : 31.0 mm			
309.0010FR	10 µL	Nat. ultra-microtip, ster.	10x 96/rack


<i>Fit Socorex®, Eppendorf®</i>			
			
Length : 45.6 mm			
302.0010FR	10 µL	Nat. microtip, sterile	10x 100/rack
302.0010AFB	10 µL	Non-sterile, autocl.	1x 1000/bag
302.0010AFR	10 µL	Non-sterile, autocl.	10x 100/rack


<i>Fit Socorex®, Eppendorf®</i>			
			
Length : 51.0 mm			
308.0020FR	20 µL	Natural tip, sterile	10x 96/rack
308.0020FS	20 µL	Single wrapped, ster.	5x 100/bag
308.0020AFB	20 µL	Non-sterile, autocl.	1x 1000/bag


<i>Fit Socorex®, Eppendorf®, Gilson®</i>			
			
Length : 69.0 mm			
308.0020GFR	20 µL	Nat. gel load. tip, ster.	4x 204/rack

<i>Fit Socorex®, Eppendorf®, Gilson®</i>			
			
Length : 54.1 mm			
308.0100FR	100 µL	Natural tip, sterile	10x 96/rack
308.0100AFB	100 µL	Non-sterile, autocl.	1x 1000/bag
308.0100AFR	100 µL	Non-sterile, autocl.	10x 96/rack

Cat. No.	Vol.	Description	Packaging
<i>Fit Socorex®, Eppendorf®, Gilson®</i>			
			
Length : 69.0 mm			
308.0100GFR	100 µL	Nat. gel load. tip, ster.	4x 204/rack

<i>Fit Socorex®, Eppendorf®, Gilson®</i>			
			
Length : 90.9 mm			
308.0200LFR	200 µL	Nat. extended tip, ster.	8x 204/rack
308.0200LFS	200 µL	Single wrapped, ster.	5x 100/bag

<i>Fit Socorex®, Eppendorf®, Gilson®</i>			
			
Length : 58.7 mm			
308.0300FR	300 µL	Natural tip, sterile	10x 96/rack
308.0300AFB	300 µL	Non-sterile, autocl.	1x 1000/bag
308.0300AFR	300 µL	Non-sterile, autocl.	10x 96/rack

<i>Fit Socorex®, Eppendorf®, Gilson®</i>			
			
Length : 88.6 mm			
309.1000FR	1000 µL	Natural tip, sterile	10x 100/rack
309.1000FS	1000 µL	Single wrapped	5x 100/bag
309.1000AFB	1000 µL	Non-sterile, autocl.	1x 1000/bag

- |   |   |                |
|---|---|----------------|
| A | = | Autoclavable   |
| B | = | Bag            |
| F | = | Filter         |
| G | = | Gel loading    |
| L | = | Long           |
| R | = | Rack           |
| S | = | Single wrapped |

For more information or for an overview of the complete Socorex product range, type

[www.socorex.com](http://www.socorex.com)

**Trade mark acknowledgement:**

Socorex®, Acura®, Calibra®, Qualitips® owned by Socorex Isba SA, Switzerland - Gilson® owned by Gilson Inc, France  
Eppendorf® owned by Eppendorf Nethler Hinz GmbH, Germany - Finnpiquette® owned by Thermo, USA

