



RLPR0516 Sparkfree Refrigerator

Part Number	RLPR0516
Capacity (Litres)	150
Style	Under Bench
Exterior Dimensions (HxWxD)	835 x 600 x 600mm
Interior Dimensions (HxWxD)	690 x 480 x 380mm
Weight	45kg
Energy Consumption	1.1kWh/24hr
Operating Temperature	+0°C to +10°C
Factory Set Point	+4°C
Interior Construction	Skinplate
Exterior Construction	Skinplate
Door Style	Self closing
Door Hinging	Right (can be reversed please ask at time of order)
Door Lock	✓
Door Open Alarm	✓
Shelves	3
Shelf dimension (WxD)	2 x 475 x 400mm 1 x 475 x 240mm
Spare shelves available	✓
Feet	4
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Min/Max Recording	✓
Access Port	✓
Off Cycle Auto Defrost	✓
Eco Friendly Refrigerant	✓
2 Year Warranty (UK only)	✓

The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this fridge.



FREE 2 YEAR WARRANTY
 LABCOLD

For extra security all sparkfree models are fitted with a built in lock and digital temperature display. The display features high and low audio and visual temperature alarms as standard in addition to a door open alarm. The controller also records the minimum and maximum temperature since the last reset.