

# Spectrophotometers at a glance



Enter this number into the search box at [www.camlab.co.uk](http://www.camlab.co.uk) to see more information, prices and live stock data

	Our Barcode	Our Barcode	Mode	Wavelength (Å)	Band Width (nm)	Data	Beam
Visible Range	1136113	Jenway 6300 Visible Spectrophotometer	Photometric Concentration	320-1000	8	PC Software Analogue RS232	Single
	1173321	Jenway 7310 Visible Spectrophotometer	Photometric Concentration Spectrum Scan Kinetics Quantitation	320-1000	5	Post scan analysis CSV,USB Pc Software Analogue RS232	Single
	1173319	Jenway 7300 Visible Spectrophotometer	Photometric Concentration	320-1000	5	Pc Software Analogue, RS232	Single
Life Science	1188507	Genova Plus Visible spectrophotometer fitted with micro-cuvette holder	Photometric Concentration Spectrum Scan Kinetics Quantitation Multi-λ	198-1000	5	Post scan analysis CSV/bmp, USB Pc Software Analogue, RS232	Single
UV / Visible Range	1136120	Jenway 6305 UV/Visible Spectrophotometer	Photometric Concentration	198-1000	8	Pc Software Analogue RS232	Single
	1173320	Jenway 7305 UV/Visible Spectrophotometer	Photometric Concentration	198-1000	5	Pc Software Analogue RS232	Single
	1173322	Jenway 7315 UV/Visible Spectrophotometer	Photometric Concentration Spectrum Scan Kinetics Quantitation	198-1000	5	Post scan analysis CSV,USB Pc Software Analogue,RS232	Single
	1189659	Jenway 6850 UV/Visible Spectrophotometer	Photometric Concentration Spectrum Scan Kinetics Quantitation Multi-λ	198-1100	0.55 1.0 2.0 4.0 5.0	Post scan analysis Pc Software USB,RS232	Double