

## Data Sheet

### UVS-26P - 6-Watt, 254nm Rechargeable UV Lamp

**Part Numbers:**

95-0187-01 (115V)                    - Non Ce Certified  
95-0187-02 (230V)                    - CE Certified

**Intensity:**

**Shortwave: 1 Tube**   - Approximately 1010  $\mu\text{W}/\text{cm}^2$  at a distance of 3" from the sample  
**Shortwave: 2 Tubes** - Approximately 1300  $\mu\text{W}/\text{cm}^2$  at a distance of 3" from the sample

**Dimensions:**

140mm Long x 76mm Wide x 249mm High

**Weight:**

1.8Kg

**Spare Parts:**

34-0013-01                                - 6-Watt, 254nm UV Tube

**Accessories:**

58-0127-02                                - 12V Adapter