

Vortex Mixers

Models

Model	Type	Control	Duty Rating	Orbit	Item No.
VXMNFS	Mini	Fixed Speed	Intermittent	4.9 mm	30392110
VXMNAL	Mini	Analog	Intermittent	4.9 mm	30392117
VXMNDG	Mini	Digital	Intermittent	4.9 mm	30392124
VXMNPS	Mini	Pulsing	Intermittent	2.5 mm	30392131
VXHDDG	Heavy-Duty	Digital	Continuous	4.9 mm	30392136
VXHDAL	Heavy-Duty	Analog	Continuous	4.9 mm	30392141
VXMPDG	Microplate	Digital	Continuous	3.5 mm	30392150
VXMPAL	Microplate	Analog	Continuous	3.5 mm	30392155
VXMTAL	Multi-Tube	Analog	Continuous	3.6 mm	30392166
VXMTDG	Multi-Tube	Digital	Continuous	3.6 mm	30392173

Accessories

Item No.	Description	VXMNFS	VXMNAL	VXMNDG	VXMNPS	VXHDDG	VXHDAL	VXMPDG	VXMPAL	VXMTAL	VXMTDG
30400151	Module Block 10 mm Test Tube	X	X								
30400152	Module Block 12/13 mm 20 Hole	X	X								
30400153	Module Block 15/16 mm 12 Hole	X	X								
30400154	Module Block 20 mm Test Tube 8 Hole	X	X								
30400155	Module Block 25 mm Test Tube 8 Hole	X	X								
30400156	Module Block Combination	X	X								
30400157	Module Block 0.5 mL Microtaper	X	X								
30400158	Module Block 6 mm Test Tube	X	X								
30400159	Module Block 1.5 mL Microtaper	X	X								
30400160	Module Block Solid	X	X								
30400161	Module Block Cuvette - 12 Cuvet	X	X								
30400162	Module Block 11.5/1.5 mL Microcentrifuge	X	X								
30400163	Module Block 250 Microliter/6 mm	X	X								
30400164	Module Block Titer Plate	X	X								
30400165	Module Block 12/13 mm 16 Holes	X	X								
30400166	Module Block Solid Double	X	X								
30400167	Module Block 4 X 35 mm Test Tube	X	X								
30400168	Module Block 50 mL Conical Centrifuge	X	X								
30400169	Module Block 0.2 mL Micro 80 Hole	X	X								
30400170	Module Block 0.2 mL Micro 64 Hole	X	X								
30400171	Module Block 0.2 mL 96 Hole Double Block	X	X								
30400172	Module Block 15 mL Conical Centrifuge	X	X								
30400173	Sand Bath, 1 Block Unit	X	X								
30400174	Sand Bath, 2 Block Unit	X	X								
30400175	Sand Bath, 4 Block Unit	X	X								
30400176	Sand Bath, 6 Block Unit	X	X								
30400177	Sand, 0.5 kg (1 Lb) Bag	X	X								
30400178	Stainless Steel Shot, 0.5 kg (1 Lb)	X	X								

Vortex Mixers

Accessories

Item No.	Description	VXMFNS	VXMNAL	VXMNDG	VXMNPS	VXHDDG	VXHDAL	VXMPDG	VXMPAL	VXMTAL	VXMTDG	
30400179	Cover Low-Temperature 2 Block	X	X									
30400180	Cover Low-Temperature 4 Block	X	X									
30400181	Cover Low-Temperature 6 Block	X	X									
30400182	Module Block For Vials 12 mm	X	X									
30400183	Module Block For Vials 15 mm	X	X									
30400184	Module Block For Vials 17 mm	X	X									
30400185	Module Block For Vials 19 mm	X	X									
30400186	Module Block For Vials 21 mm	X	X									
30400187	Module Block For Vials 23 mm	X	X									
30400188	Module Block For Vials 25 mm	X	X									
30400189	Module Block For Vials 28 mm	X	X									
30400190	Module Block For Vials 16 mm	X	X									
30400191	Module Block, 2.0 mL Tube	X	X									
30400192	Module Block, 2.0 mL Corning Tube	X	X									
30400193	Module Block Combo, 50 mL, 15 mL, 1.5 mL	X	X									
30400194	Module Block Combo - Micro Tube	X	X									
30400195	Module Block 12 Holes 17-18 mm	X	X									
30400246	Stainless Steel Probe Kit	X	X									
30400196	8 cm Head with Rubber Cover	X	X	X	X							
30400197	Tube Holder Kit with Retainer		X	X	X							
30400198	Micro Tube and Well Plate Holder Kit		X	X	X							
30400199	Micro Well Plate Kit with Retainer		X	X	X							
30400200	Microtube Kit with Retainer		X	X	X							
30400201	Foam Insert Variety Pack		X	X	X							
30400202	Single Tube Holder		X	X	X							
30400203	Single Tube Holder Adapter Vortex Genie		X	X	X							
30400204	0.5 mL Micro-Tube Holder		X	X	X							
30400205	1.5-2.0 mL Micro-Tube Holder		X	X	X							
30400206	1.5-2.0 mL Stainless Steel Tube Holder		X	X	X							
30400207	Ampule Tube Holder - Small (15-17mm)		X	X	X							
30400208	Ampule Tube Holder - Large (10-17mm)		X	X	X							
30400209	Foam Insert Blank					X	X					
30400210	Cup Head, VXHD and VXMP Vortex Mixers					X	X	X	X			
30400211	Large Vessel Holder For 500 ml/1 L Flask					X	X					
30400212	4 Place Stackable Microplate Holder					X	X					
30400213	Double Microplate Holder					X	X	X	X			
30400214	4 Place Microplate Holder					X	X					
30400215	Single Microplate Holder					X	X	X	X			
30400216	Foam Insert For (52) 0.5 mL Microtubes					X	X					
30400217	Foam Insert (38) 1.5-2.0 mL Microtubes					X	X					
30400218	Small Vessel Holder 125 & 250 mL Flasks					X	X					
30400219	Single Tube Holder, VXHD Vortex Mixers					X	X					
30400220	Foam Insert For (34) 12-13 mm Tubes					X	X					
30400221	Foam Insert For (20) 15-18 mm Tubes					X	X					
30400222	Foam Insert For (18) 19-21 mm Tubes					X	X					
30400223	Foam Insert For (13) 22-25 mm Tubes					X	X					
30400224	Foam Insert For (4) 26-29 mm Tubes					X	X					

Vortex Mixers

Accessories

Item No.	Description	VXMFNS	VXMNAL	VXMNDG	VXMNPS	VXHDDG	VXHDAL	VXMPDG	VXMPAL	VXMTAL	VXMTDG	
30400225	Cover For Universal Holder					X	X					
30400226	Universal Holder					X	X					
30400238	10 mm Tube Rack (Gray)									X	X	
30400239	12 mm Tube Rack (Blue)									X	X	
30400240	13 mm Tube Rack (Yellow)									X	X	
30400241	16 mm Tube Rack (Green)									X	X	
30400242	50 mL Tube Rack (Red)									X	X	
30400243	25 mm Tube Rack (White)									X	X	
30400244	Post Extension Kit									X	X	
30400245	Tray Pad Set - Top And Bottom									X	X	
30400227	Insert Retainer		X	X	X							
30400228	Vessel Harness, Pack of 2		X	X	X							
30400229	Holder For 9-13 mm Tubes, Pack of 2		X	X	X							
30400230	Holder For 14-19 mm Tubes, Pack of 2		X	X	X							
30400231	Holder For 20-25 mm Tubes, Pack of 2		X	X	X							
30400232	Micro-Tube Holder, Pack of 2		X	X	X							
30400233	Micro-Well Plate Holder, Pack of 2		X	X	X							
30400234	Flat Foam Insert, Pack of 2		X	X	X							
30400235	Cup Head		X	X	X							
30400236	7.6 cm Head (3")	X	X	X	X							
30400237	7.6 cm Rubber Head Cover (3")	X	X	X	X							