

C2006 Fixed angle up to 6000RPM

“functional simplicity”

- Speed** 500-6000Rpm
in 10 Rpm steps
- Max Rcf** 3800G
- Timer** 0-99Minutes
and constant
- Dims HWD** 27x33x46Cms
- Weight** 16.5Kilos
without rotor



Rotors Available



BR5412

max speed 6000 Rpm 3,800G
12 x 15ml tube, fixed angle rotor
Tube size maximum H120 x W17 mm



RM10

ADAPTORS for BR5412 (above) 5-10ml
Tube size maximum H100 x W13 mm



BR5406

max speed 6000 Rpm 3,800G
6 x 50ml Tube, fixed angle rotor
Tube size maximum H115 x W30mm



RM25

ADAPTORS for BR5406 (above) 25ml
Tube size maximum H100 x W25 mm

RM15

ADAPTORS for BR5406 (above) 15ml
Tube size maximum H120 x W17 mm