

EBA 270 TABLETOP CENTRIFUGE

FOR **CLINICAL** APPLICATIONS

FEATURES:


- For GEL BLOOD TUBES
- Ideal for Physicians Offices
- 90° Swing-out Built-in Rotor
- Maintenance free, brushless motor
- Run-Silent™ Operation (≤ 60 dB, softer than normal conversation)
- Safety-Lock™ Lid
- Five (5) Year Warranty
- FDA registered centrifuge manufacturer

Hettich
CENTRIFUGES



SMALL IN SIZE, BIG ON RESULTS

Blood separation on the EBA 270 spinning a standard 15 mL Gel Tube.



(Built in 90° horizontal rotor)



NEW

EBA 270 **BLOOD TUBE** PACKAGE



Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
6	1.6 - 7	13 x 75/ 100	4,000	1,807/ 2,254	(1) Centrifuge (1) Swing-out Rotor (6) Adapters for 13 x 75 mm (6) Adapters for 13 x 100 mm
	8.5 - 10	16 x 75/ 100			
Pkg No.:			270BLOOD		

DISTRIBUTED BY:

Distributed and Supported in the UK by:



24 Norman Way Industrial Estate, Over, Cambridge, CB24 5WE
Tel: +44(0)1954 233 100 Fax: +44(0)1954 233 101
Email: sales@camlab.co.uk Web: www.camlab.co.uk

DIMENSIONS (HxWxD) 9.4" x 12.8" x 15.31"
TUBE TYPES ROUND BOTTOM, BLOOD

MAX. CAPACITY:
6 x 15 mL

Copyright © 2012 Hettich Instruments, LP. All rights reserved. For the most up to date information on Hettich please visit hettweb.com. Printed in the USA using soy based inks on FSC certified recycled paper (25% pre-consumer waste & 25% post consumer waste). PS-EBA-270-E RevA 01-2013