



Gay Lussac - Alcoholometer

250mm long, measuring alcohol % by volume and adjusted at 15°C (60°F).

Product Specifications:

Material:	Soda glass, steel shot ballast, Synthetic wax resin
Range:	0 to 100
Accuracy:	+/- 2
Divisions:	1
Bulb Diameter	16mm
Overall Length:	250mm
Weight:	29g (including protective plastic tube)



Description	Range % Alcohol	Divisions	Product No
Alcoholometer	0 to 100	1	10/670/0

PDS 0518:13:R0