



KLEENGUARD® Product Description

KLEENGUARD® A40 Overshoe with Sole

- Code: 98710 (M & L) and 98720 (XL & XXL) (200 items / case)
- Description: An overshoe with elasticated ankle and PVC sole
White
- Protection: Although this item of equipment is manufactured from a fabric which has liquid splash and dust barrier properties (suitable for Type 5/6), it is not intended as primary protection to the foot, and therefore is marked CE simple.

Approximate Dimensions

	Size	M&L	XL & XXL
Foot:	A ⇄ B	33.5 cm	38 cm
Elastic:	C ⇄ D	12.5 cm	12.5 cm
Height:	G ⇄ H	13.3 cm	13.3 cm

