













### BD Vacutainer® EDTA Tubes

- Single use, evacuated, sterile blood collection tube
- Contains K2EDTA or K3EDTA
- Used to obtain a whole blood or EDTA plasma sample
- Pack/Case Quantity: 100/1000
- Centrifugation: N/A

	<i>Cat. No.</i>	<i>Draw (ml)</i>	<i>Size (mm)</i>	<i>Additive</i>	<i>Material</i>	<i>Label</i>	<i>Closure</i>	<i>Cap Colour</i>
	368841	2	13x75	K2EDTA	PET	Paper	Hemogard™	Translucent Lavender
	368274	2	13x75	K2EDTA	PET	See Through	Hemogard™	Translucent Lavender
	368856	3	13x75	K2EDTA	PET	Paper	Hemogard™	Translucent Lavender
	368499	3	13x75	K2EDTA	PET	See Through	Hemogard™	Translucent Lavender
	368861	4	13x75	K2EDTA	PET	Paper	Hemogard™	Lavender
	367839	4	13x75	K2EDTA	PET	Block	Hemogard™	Lavender
	367862	4	13x75	K2EDTA	PET	See Through	Hemogard™	Lavender
	367864	6	13X100	K2EDTA	PET	Paper	Hemogard™	Lavender
	365900	6	13X100	K2EDTA	PET	See Through	Hemogard™	Lavender
	367525	10	16X100	K2EDTA	PET	Paper	Hemogard™	Lavender
	367836	2	13x75	K3EDTA	PET	Block	Hemogard™	Translucent Lavender
	367858	2	13x75	K3EDTA	PET	See Through	Hemogard™	Translucent Lavender
	368857	3	13x75	K3EDTA	PET	Block	Hemogard™	Translucent Lavender
	368860	4	13x75	K3EDTA	PET	Block	Hemogard™	Lavender
	368270	4	13x75	K3EDTA	PET	See Through	Hemogard™	Lavender