

# Brannan Product Data Sheet

PDS0511.20.REV1

## SPH004 - Special Purpose Specific Gravity Hydrometer - Specific Gravity & Quevenne Lactometers

These Lactometers are 160mm long milk hydrometers, adjusted at 15.6°C (60°F). Scales calibrated in specific gravity or Quevenne. 1° Quevenne equals 0.001°SG. All hydrometers can be supplied with UKAS certification, certificate traceable to Recognised National Standard or Statement of Conformity, if required.



### Product Features:

- Adjusted at 15.6°C (60°F)
- 160mm long
- Made in UK
- Supplied in protective, resealable plastic tube

### Product Specifications:

<b>Material:</b>	Soda glass, steel shot ballast, synthetic wax resin
<b>Ranges:</b>	See table below
<b>Accuracy:</b>	See table below
<b>Divisions:</b>	See table below
<b>Overall Length:</b>	160mm
<b>Bulb Diameter:</b>	16mm
<b>Approx Weight:</b>	22g (including 8g protective plastic tube)

Range	Accuracy at 20°C	Division	Product Number
1.025 to 1.035°SG	+/- 0.002	0.001	10/601/0
15 to 40° Quevenne/ 1.015 to 1.025°SG	+/- 2.000	1.000/ 0.100	10/660/0

[www.brannan.co.uk](http://www.brannan.co.uk)

**Measurement for Life.**

