



Hydrometer Cylinders

Soda glass hydrometer cylinders available with a standard internal diameter of 37mm and plastic foot.



Internal Depth	Suitable for Hydrometer Series	Product No
200mm	S50	10/000/0
300mm	M50 & M100	10/000/1
330mm	L50	10/000/2

PDS 0218:12:R0

Please note: all information on this data sheet is correct at the time of printing, but our policy is one of continuous product development. We reserve the right to change specification without notice.

Brannan Thermometers & Instrumentation

Leconfield Industrial Estate • Cleator Moor • Cumbria • CA25 5QE • UK

tel: +44 (0)1946 816600 fax: +44 (0)1946 816625 email: sales@brannan.co.uk web: www.brannan.co.uk



0582