



Camlab OFFERS

14th October - 14th December

Hach DR2800
New 16mm
version now
available!

EXCLUSIVE OFFER!
ONLY £1750
including **FREE Service Voucher** and a **Socorex pipette worth £113**

~~RRP £2,550~~
YOU SAVE
£800



The Hach DR2800 water testing spectrophotometer is now available from Camlab in a new 16mm format which is compatible with most manufacturers test systems – no need to change trusted chemistries.

Barcode recognition is available on the 13mm cuvette tests – no mistakes as the test parameters are automatically recognised.

DR2800 for 16mm vials, on-site training, FREE service voucher and a FREE reagent sample pack

www.camlab.co.uk/offers/en677

As above but for the 13mm vial system

www.camlab.co.uk/offers/en678

+ FREE REAGENTS SAMPLE PACK

More Offers Inside

up to 50% off

Prize Draw - Win an Apple iPod 30 Gb video



Camlab Limited, Camlab House, Norman Way, Over, Cambridge, CB4 5WE
Tel: 01954 233110 • Fax: 01954 233101
email: sales@camlab.co.uk • www.camlab.co.uk

up to
**50%
OFF**

Camlab CUB Glass Washers

These fabulous offers are only available as a limited offer to launch these two new stainless steel washers from Camlab. The **CUB Compact** is a slimline unit for use by the lower through-put lab who still require high performance. The **CUB Easy** is robust yet still fits under the bench for the larger lab or core facility who demand excellence with simplicity.

**FREE
DELIVERY
ON ALL
OFFERS**

Both units are quiet and can be installed in the lab and the pre programming makes them easy to use. They operate on 1 or 2 levels and are suitable for a wide variety of glassware. All offers include 40L of neutralising liquid and commissioning which covers checking the installation and operation as well as user training.



CUB Compact
450x600x850mm

CUB Compact + 2 std racks RRP ~~£4144~~ now only £2240

www.camlab.co.uk/offers/en673

CUB Compact + Jet Rack RRP ~~£4536~~ now only £2685

www.camlab.co.uk/offers/en671

**+ FREE 40L NEUTRALISING LIQUID
+ FREE COMMISSIONING**

CUB Easy + 2 std racks RRP ~~£5053~~ now only £2985

www.camlab.co.uk/offer/en674

CUB Easy + Jet Rack RRP ~~£5696~~ now only £3195

www.camlab.co.uk/offer/en672




CUB Easy 600x600x850mm

The standard racks hold a variety of beakers and flasks on wire pegs. For narrow necked flasks the Jet Rack has specific supports with injection nozzles for injecting water and detergent directly into the glassware.



Std Rack 1



Jet Rack



Std Rack 2

Camlab COD Colorimeter offer

Camlab is able to offer you COD vials in a choice of concentrations and vial sizes for immediate delivery from our UK stock. Whether you choose Hach or Camlab's own brand reagents you can be sure of rapid delivery. Place your next order with us before 4th December 2006 to get a free CW2020 colorimeter from the new Camlab Water range.

Order any 16 packs of COD vials supplied by Camlab and get a free Camlab Water CW2020 colorimeter.

www.camlab.co.uk/offers/en675



**FREE CW2020
COLORIMETER
WORTH £345**

Hettich Microcentrifuge

Small and extremely powerful, the Hettich Mikro 120 is compact and provides quiet and dependable centrifugation of micro tubes up to 14,000 rpm (17.553 xg). As it has an easy to use touch panel display, a 1 to 99 min timer, an impulse button and automatic lid lock.

Mikro 120, 24 x 1.5/2.0ml place rotor, 1000 x 1.5ml tubes, 3 year warranty £720

www.camlab.co.uk/offers/en679

- **FREE ROTOR WORTH £160**
- **1000 FREE TUBES WORTH £27**
- **3 YEAR WARRANTY WORTH £54**



Package only £720

RRP ~~£960~~



Camlab COD Vial Disposal offer

Camlab is able to offer you many types and brands of COD vial, but realise disposal is also an important consideration. To help you stretch your budget towards 2007 we have put together several offers to help maximise your testing.

Buy 12 packs of any COD vial supplied by Camlab and receive a free COD vial disposal (25 litre drum) www.camlab.co.uk/offers/en676


How the Camlab COD Disposal works - We send you a safety disposal drum. When you have filled your drum with vials simply call us and we will arrange for collection by authorised carrier and carry out proper disposal of your waste.



**FREE COD DISPOSAL
SAVE £170**



Ensure that your COD vials are disposed of safely and in accordance with the Duty of Care section of the Environmental Protection Act (1990) by using the Camlab Vial Disposal Service. Each producer of hazardous waste has to register with the Environment Agency. Register at www.environment-agency.gov.uk or by phoning 08708 502858



Order online at www.camlab.co.uk

FAX to 01954 233 101

To get your FREE Camlab pen set and your chance to win an IPOD 30GB video simply complete the quick survey below and fax back to us.

Who is your usual supplier for the following...

- a. Water Testing Equipment. _____
 b. Water Testing Reagents. _____
 c. Laboratory Equipment. _____
 d. Laboratory Reagents. _____
 e. Repairs and Maintenance. _____
 of the above equipment.

How would you rate Camlab for supply of these items compared to your usual supplier?

- Much Better Better Same Worse Don't Know
 Much Better Better Same Worse Don't Know
 Much Better Better Same Worse Don't Know
 Much Better Better Same Worse Don't Know
 Much Better Better Same Worse Don't Know



Does your company purchase the following support services for laboratory or field testing equipment?

- a. Service contracts
 Yes No Don't Know
 b. Extended warranties
 Yes No Don't Know

If you answered 'No' do you think your company is likely to purchase these in the next 24 months?

- a. Service contracts
 Yes No
 b. Extended warranties
 Yes No

RAPID REPLY Fax Form – 01954 233 101 or Phone 01954 233110

Title: _____ First name: _____
 Surname: _____
 Department: _____
 Company: _____
 Division: _____
 Address: _____

 City/Town: _____
 County: _____ Postcode: _____
 Phone: _____
 Fax: _____
 Email: _____

Data protection: Camlab may contact you from time to time with news and offers. If you would prefer not to receive any email communications from us please tick here
 If you would prefer not to receive any subsequent posted information from us please tick here If you would like your details removed from our mailing list entirely please tick here Thank you.

TICK HERE

- Please send me my FREE pen pack
 Please enter me in the draw for an IPOD 30GB
 Please send me a FREE CUB Warranty voucher
 Please send me FREE COD vial sample
 Please send me more information on...
 Hach DR2800 CUB Washers COD Vials
 COD Vial Disposal Service Hettich Microcentrifuge
 Please send me a Camlab catalogue
 Environment Laboratory
 Please keep me up to date with future offers by email



EXPRESS ORDER FORM: Please process my order as detailed below ...

Offer No.	Description	Offer Price £	QTY Required	Total price - £
EN 671	CUB Compact + Jet Rack + Commissioning + Neutraliser	£2,685.00		
EN 672	CUB Easy + Jet Rack + Commissioning + Neutraliser	£3,195.00		
EN 673	CUB Compact + STD Racks + Commissioning + Neutraliser	£2,240.00		
EN 674	CUB Easy + STD Racks + Commissioning + Neutraliser	£2,985.00		
EN 677	Hach DR2800 + FREE Soccorex pipette (16mm System)	£1,750.00		
EN 678	Hach DR2800 + FREE Soccorex pipette (13mm System)	£1,750.00		
EN 679	Hettich Mikro 120 + FREE Rotor & 1000 1.5ml Tubes	£720.00		
EN 675	16 x pack of COD Vials + FREE CW2020 colorimeter	£		
EN 676	12 x pack of COD Vials + Free COD Vial disposal Service	£		

Order Number	
Camlab Account Number	

SubTotal	
Delivery	FREE
VAT @ 17.5%	
Grand Total	