

Camlab Plastics

Blackline Light-safe Storage Solutions



0.5ml & 1.5ml

5ml

15ml

50ml

A RANGE OF STORAGE SOLUTIONS FOR LIGHT SENSITIVE SAMPLES

Camlab offers a new family of products to keep your light sensitive samples, reagents and solutions safe. Unlike amber tinted tubes, these products are completely opaque therefore protecting even the most light sensitive compounds. Camlab Blackline tubes and boxes eliminate the inconvenience of having to wrap the tubes in aluminium foil to protect samples. Ideal for use with fluorescent probes or ethidium bromide.

The entire range is autoclavable and features storage boxes, 0.5 and 1.5ml centrifuge tubes, our unique 5ml preparaton tube as well as standard 15 and 50ml conical centrifuge tubes.

For further details, samples and pricing please contact us using the contact information overleaf.



Features





- NO MORE TINFOIL!
- For the protection of your photosensitive samples and reagents
- 0.5ml, 1.5ml and 5ml tubes are ideal for small volume samples
- 100 place box for the storage of small volume tubes
- Extended range now includes 15ml and 50ml tubes
- All products from the range are autoclavable
- High RCF rating on all tubes



camlab

www.camlab.co.uk

ORDERING INFORMATION

Catalogue no.	Description	Pack Qty	Case Qty
	0.5ml & 1.5ml microcentrifuge tubes Made from medical grade polypropylene, Blackline micro centrifuge tubes can be stored at temperatures down to -100°C and centrifuged at up to 15,000xg. Click-fit cap provides a tight seal to prevent opening during centrifugation. Autoclavable.		
RTP/77099-BK	0.5ml Blackline microtube, PP	500	10000
RTP/77120-BK	1.5ml Blackline microtube, PP	500	10000
RTP/77110-BK	1.5ml Blackline microtube, PP	1000	10000
RTP/77100-BK	1.5ml Blackline microtube, PP	5000	25000
	100 Position Box These stackable 100 place boxes with hinged lids are suitable for use with 0.5, 1.5 & 2.0ml microtubes or cryogenic vials up to 2ml. The lids have an alpha numeric grid for easy sample location and matt surfaces on the lid to allow easy attachment of labels. Autoclavable.		
RTP/72101-BK	Blackline 100 place storage box.	1	100
	5ml Conical preparation tube. Made from medical grade polypropylene, Blackline preparation tubes are ideal as an intermediate sized tube and can be stored at temperatures down to -100°C and centrifuged up to 5,000xg. Click-fit cap provides tight seal to prevent opening during centrifugation. Autoclavable.		
RTP/0310-BK	5ml Blackline centrifuge tube, PP	250	5000
	15 & 50ml Centrifuge tubes. Made from polypropylene with a leakproof screwcap. These robust tubes can be centrifuged up to 8,400xg and 10,000xg respectively. They also have a frosted writing area for easy permanent identification. Autoclavable.		
RTP/1500 -BK	15ml Blackline centrifuge tube, PP	50	500
RTP/5000-BK	50ml Blackline centrifuge tube, PP	50	500

ALSO AVAILABLE FROM CAMLAB




Double Tamper Evident Bottle



Micryo Cryogenic Labels




100 Position Freezer Boxes



Tamper Evident 1.5ml Tube



1.5ml Graduated Micro Tube



Camlab Limited, Camlab House, Norman Way, Over, Cambridge, CB24 5WE
 Tel: +44 (0)1954 233110 • Fax: +44 (0)1954 233101
 email: plastics@camlab.co.uk • www.camlab.co.uk