

Certificate of Analysis List

For request number 419307

Catalog Number Entered	Lot Number Entered	Related Catalog Number	Related Lot Code	Description
2659949	1223	N/A	N/A	StablCal [®] Standard, 10 NTU

Total Enclosures: 1



An ISO 9001 Certified Company

Certificate of Analysis

Page 1

COMMODITY: **StablCal[®] Standard, 10 NTU**COMMODITY NUMBER: **26599-49**

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: **A1223****8/11/2011****8/15/2011**

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	9.5 to 10.5 NTU	10.23 NTU

The expiration date is Aug 2013

Formazin and StablCal[®] solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal[®] as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

A handwritten signature in cursive script that reads "Scott Als".

Certified by _____

Scott Als
Analytical Services Chemist