

DURAN® RANGE OF GLS 80® BOTTLES AND CAPS






Distributed By: Camlab Ltd
Unit 24, Norman Way Industrial Estate
Over, Cambridge, CB24 5WE, United Kingdom
T: +44 (0) 1954 233 110 E: sales@camlab.co.uk






**DURAN
WHEATON
KIMBLE**

Excellence in your hands

Product range		DURAN®GLS 80®WIDE MOUTH LABORATORY BOTTLE		DURAN®GLS 80®WIDE MOUTH LABORATORY BOTTLE, AMBER		DURAN®GLS 80®WIDE MOUTH LABORATORY BOTTLE, PROTECT			
Borosilicate 3.3. glass bottle body									
Caps ¹		Available with or without GLS80®screw cap from PP		Available with or without GLS80®screw cap from PP		Available with or without GLS80®screw cap from PP			
Temperature resistance		Bottle: -70°Cto +500°C Cap: -40°Cto +140°C		Bottle: -70°Cto +500°C Cap: -40°Cto +140°C		Bottle: -30°Cto +135°C Cap: -40°Cto +140°C			
Main advantage		<ul style="list-style-type: none"> • Tried and tested, DURAN®glass • Wide mouth enables easy access 		<ul style="list-style-type: none"> • USP <660> and USP <671> (Spectral Transmission) compliant • UVprotection up to ca. 500nm 		<ul style="list-style-type: none"> • The coating provides scratch, leak and splinter protection² • UVprotection up to ca. 380nm (ca. 500nm Amber) 			
Colour		clear		amber		clear		amber	
mL	GLS thread	With screw cap (PP)	Without screw cap	With screw cap (PP)	Without screw cap	With screw cap (PP)	Without screw cap	With screw cap (PP)	Without screw cap
250	80	21 860 36 56	21 860 36 07	21 866 36 53	21 866 36 04	21 865 36 52	21 865 36 03	-	-
500	80	11 126 27	11 783 92	11 601 46	11 784 29	11 601 52	-	11 673 08	21 866 44 36
1000	80	11 127 13	11 784 24	11 601 47	11 784 30	11 601 63	-	11 673 09	21 866 54 32
2000	80	11 127 15	11 784 25	11 601 48	11 784 31	11 601 64	-	-	21 866 63 34
3500	80	21 860 69 58	21 860 69 04	21 866 69 59	21 866 69 01	21 865 69 53	21 865 69 09	-	-
5000	80	11 139 49	11 784 26	11 601 49	11 784 32	11 601 65	-	-	21 866 73 39
10000	80	11 139 50	11 784 27	11 601 50	11 784 33	-	-	-	-
20000	80	11 139 51	11 784 28	11 601 51	11 784 34	-	-	-	-

¹ All these bottles are compatible with the full range of DURAN®GLScaps, including chemically resistant, venting membrane, pharmaceutical grade, and connection system caps.

² For the bottle sizes from 250 – 5 000 mL

DURAN®GLS 80®QUICK RELEASE CLOSURE FROM PP		DURAN®GLS 80®QUICK RELEASE CLOSURE WITH CAP LINER (PSU COMPOUND)			DURAN®GLS 80®MEMBRANE VENTED SCREW CAP FROM PP
					
A matching PPs pouring ring is available		A matching PTFE pouring ring is available			Use with either PP or PTFE pouring ring
-40 °C to +140 °C		-45 °C to +180 °C			-40 °C to +140 °C
Permits opening and closing of the DURAN® GLS80® bottle with only a three-quarter turn		The PSU material with PTFE coated liners offers improved chemical, thermal and mechanical properties			Ideal for autoclave application, membrane permits pressure equalisation
blue		white			blue
Screw cap	Pouring ring	Screw cap	Pouring ring	Replacement cap liner	Screw cap
11 127 16	11 601 66	11 658 88	11 673 07	11 529 21	29 118 91 05