

Distributed By: Camlab Ltd
Unit 24, Norman Way Industrial Estate
Over, Cambridge, CB24 5WE, United Kingdom
T: +44 (0) 1954 233 110 E: sales@camlab.co.uk



DURAN
WHEATON
KIMBLE

Excellence in your hands

DURAN®PURE

DEVELOPED FOR PHARMA

DURAN®PURE is a range of high-quality bottles and closures developed for the needs of the pharmaceutical industry. DURAN®PURE meets the requirements of the international pharmacopeias and is available with the certificates required by the regulatory authorities.

FEATURES

- EP/USP Type I borosilicate glass
- Size range 25 ml - 20 l
- Protective cover prevents contamination
- Single batch delivery CoA available
- 100% traceability
- 100% visual inspection
- Change control
- Drug Master File

APPLICATION AREAS

- Manufacturing of pharmaceutical products
- Storage of bulk drug intermediates
- Transportation of APIs



DURAN®PUREBottles

| Article number clear glass | Article number amber coloured | Volume [ml] | DIN thread GL | Bottles per box [pcs.] | Bottles per pallet [pcs.] |
|-------------------------------|----------------------------------|----------------|------------------|---------------------------|------------------------------|
| 81 801 14 01 | 81 806 14 06 | 25 | 25 | 10 | 940 |
| 81 801 24 06 | 81 806 24 02 | 100 | 45 | 10 | 540 |
| 81 801 36 08 | 81 806 36 04 | 250 | 45 | 10 | 280 |
| 81 801 44 07 | 81 806 44 03 | 500 | 45 | 10 | 360 |
| 81 801 54 03 | 81 806 54 08 | 1 000 | 45 | 10 | 210 |
| 81 801 63 05 | 81 806 63 01 | 2 000 | 45 | 10 | 120 |
| 81 801 73 01 | 81 806 73 06 | 5 000 | 45 | 1 | 48 |
| 81 801 86 06 | 81 806 86 02 | 10 000 | 45 | 1 | 30 |
| 81 801 91 05 | 81 806 91 01 | 20 000 | 45 | 1 | 12 |

ALSO AVAILABLE: DURAN®PUREGLS80®bottles and screw caps

DURAN®PUREPremium Screw Caps

DURAN®PUREPremium Screw Caps made from PFA (fluoropolymer) are chemically and thermally very resistant and compatible with autoclaving and other common sterilization methods. The products include silicone seals which are coated with highly inert PTFE.

FEATURES

- Caps made from PFA with PTFE coated silicone seal
- Produced and packaged in ISO Class 7 clean rooms
- Single batch delivery CoA available
- 100% traceability
- Change control
- Drug Master File



DURAN®PUREPremium Screw Cap



DURAN®PUREPremium Lip Seal Cap

| Article number | DIN thread GL | Description | Screw Caps per PE bag [pcs.] |
|----------------|------------------|--------------------------------|---------------------------------|
| 89 108 13 06 | 25 | DURAN®PUREPremium Screw Cap | 5 |
| 89 108 86 79 | 45 | DURAN®PUREPremium Screw Cap | 5 |
| 89 886 28 02 | 45 | DURAN®PUREPremium Lip Seal Cap | 5 |



DWK Life Sciences GmbH
Hattenbergstraße 10
55122 Mainz, Germany
www.DWK-LifeSciences.com

