

Distributed By: Camlab Ltd
Unit 24, Norman Way Industrial Estate
Over, Cambridge, CB24 5WE, United Kingdom
T: +44 (0) 1954 233 110 E: sales@camlab.co.uk



**DURAN
WHEATON
KIMBLE**

Excellence in your hands

DURAN®PURE GLS 80®

DEVELOPED FOR PHARMA

DURAN®PURE is a range of high-quality bottles and closures developed for the needs of the pharmaceutical industry. DURAN®PURE meets the requirements of the international pharmacopeias and is available with the certificates required by the regulatory authorities.

FEATURES

- EP/USP Type I borosilicate glass
- Size range 500 ml - 20 l
- GLS80® thread for fast closing & opening
- Protective cover prevents contamination
- Single batch delivery CoA available
- 100% traceability
- 100% visual inspection
- Change control
- Drug Master File

APPLICATION AREAS

- Manufacturing of pharmaceutical products
- Storage of bulk drug intermediates
- Transportation of APIs



DURAN®PURE GLS 80®Bottles

| Article number clear glass | Article number amber coloured | Volume [ml] | Bottles per box [pcs.] | Bottles per pallet [pcs.] |
|-------------------------------|----------------------------------|----------------|---------------------------|------------------------------|
| 81 860 44 03 | 81 866 44 09 | 500 | 10 | 280 |
| 81 860 54 08 | 81 866 54 05 | 1 000 | 10 | 210 |
| 81 860 63 01 | 81 866 63 07 | 2 000 | 10 | 120 |
| 81 860 73 06 | 81 866 73 03 | 5 000 | 1 | 48 |
| 81 860 86 02 | 81 866 86 08 | 10 000 | 1 | 30 |
| 81 860 91 01 | 81 866 91 07 | 20 000 | 1 | 12 |

ALSO AVAILABLE: DURAN®PURE bottles and screw caps with GL45 thread

DURAN®PURE GLS 80®Screw Cap

The DURAN®PURE GLS 80® Screw Cap made from PSU is chemically and thermally very resistant and compatible with autoclaving and other common sterilization methods. The product includes a silicone seal which is coated with highly inert PTFE.

FEATURES

- Cap made from PSU with PTFE coated silicone seal
- **GLS 80® thread for fast closing & opening**
- Single batch delivery CoA available
- 100% traceability
- Change control



| Article number | Screw Caps per PE bag [pcs.] | Screw Caps per pallet [pcs.] |
|----------------|---------------------------------|---------------------------------|
| 89 116 91 09 | 5 | 1435 |



DWK Life Sciences GmbH
Hattenbergstraße 10
55122 Mainz, Germany
www.DWK-LifeSciences.com

