

RANGE OF BULB PIPETTES FROM SODA-LIME GLASS



Product Range	BULB PIPETTES		
Accuracy class	Class AS		Class B
Material of the pipettes: soda-lime glass (AR®glass)			
Certificate*	Batch certificate	Batch certificate	–
Maximum recommended temperature without affecting accuracy (at 20 °C)	250 °C	250 °C	250 °C
Print color	Amber stain graduation	Blue	Amber stain graduation
ml			
0.5	–	23 339 00 51	–
1	24 338 01 09	23 339 01 05	24 337 01 08
2	24 338 02 03	23 339 02 08	24 337 02 02
3	–	23 339 03 02	–
4	–	23 339 04 05	–
5	24 338 07 09	23 339 05 08	24 337 07 08
6	–	23 339 06 02	–
7	–	23 339 07 05	–
8	–	23 339 08 08	–
9	–	23 339 09 02	–
10	24 338 08 03	23 339 10 07	24 337 08 02
15	–	23 339 15 04	–
20	24 338 12 08	23 339 20 03	24 337 12 07
25	24 338 14 05	23 339 25 09	24 337 14 04
30	–	23 339 30 08	–
40	–	23 339 40 04	–
50	24 338 17 05	23 339 50 09	24 337 17 04
100	24 338 24 01	23 339 00 02	24 337 24 09

* Batch certificates also available online at www.DWK-LifeScience.com

AR®glass = registered trademark of SCHOTT AG

DWKLife Sciences GmbH • Hattenbergstraße 10 • 55122 Mainz • Germany • www.DWK-LifeSciences.com

Phone: +49 6131-1445 4131 • Fax: +49 6131-1445 4016 • sales@DWK-LifeSciences.com

RANGE OF MEASURING PIPETTES FROM SODA-LIME GLASS



Product Range	MEASURING PIPETTES						
Accuracy class	Class AS					Class B	
Material of the pipettes: soda-lime glass (AR®glass)							
Certificate*	Batch certificate	Batch certificate	Batch certificate	Batch certificate	Batch certificate	—	—
Maximum recommended temperature without affecting accuracy (at 20 °C)	250 °C	250 °C	250 °C	250 °C	250 °C	250 °C	250 °C
Print color	Amber stain graduation	Blue	Amber stain graduation	Blue	Blue	Amber stain graduation	Amber stain graduation
ml	TYPE 3	TYPE 1	TYPE 2	TYPE 2	TYPE 3	TYPE 1	TYPE 3
0.1	—	—	—	—	—	24 343 01 02	24 344 01 03
0.2	—	—	—	—	—	24 343 03 08	24 344 03 09
0.5	—	23 346 06 06	23 347 06 07	23 348 06 08	23 349 06 09	24 343 06 08	24 344 06 09
1	24 345 11 09	23 346 11 05	23 347 11 06	23 348 11 07	23 349 11 08	24 343 11 07	24 344 11 08
2	24 345 17 09	23 346 16 02	23 347 16 03	23 348 16 04	23 349 16 05	24 343 16 04	24 344 16 05
5	24 345 23 02	23 346 23 07	23 347 23 08	23 348 23 09	23 349 23 01	24 343 23 09	24 344 23 01
10	24 345 29 02	23 346 29 07	23 347 29 08	23 348 29 09	23 349 29 01	24 343 29 09	24 344 29 01
20	—	23 346 32 09	23 347 32 01	23 348 32 02	23 349 32 03	—	—
25	24 345 34 01	23 346 34 06	23 347 34 07	23 348 34 08	23 349 34 09	24 343 34 08	24 344 34 09
50	—	23 346 36 03	23 347 36 04	23 348 36 05	23 349 36 06	—	—

Distributed By: Camlab Ltd
Unit 24, Norman Way Industrial Estate
Over, Cambridge, CB24 5WE, United Kingdom
T: +44 (0) 1954 233 110 E: sales@camlab.co.uk

TYPE1— partial delivery, zero point at the top
TYPE2— total delivery, nominal volume at the top
TYPE3— total delivery, zero point at the top

