

Velocity 6 μ

Micro Fuge



Features & Characteristics

- > Microtubes rotor and 8-strip rotor included
- > Rotors can be easily exchanged without tools
- > Compact design
- > Simple & user-friendly operation
- > Quiet operation (noise level \leq 45dB)

VELOCITY 6 μ SPECIFICATIONS	Velocity 6 μ
Maximum Speed	6,000 RPM
Maximum RCF	2,910 x g
Run Time	Continuous
Rotor (Standard)	8 x 1.5/2.0 ml, 2 x 8-strip PCR or 16 x 0.2 ml PCR tubes
Rotor (Optional)	-
Dimensions	155 x 175 x 125 mm (w x d x h)
Weight	1 kg
Power Supply	AC 110-220V 50/60Hz

Velocity 6 μ 6 μ Ordering Information

PRODUCT	CATALOG NUMBER#
Velocity 6 μ Micro Fuge, with microtube rotor & tube-stripes rotor, 110-220V, 50/60Hz	V6U

Distributed and Supported in the UK by:



24 Norman Way Industrial Estate, Over, Cambridge, CB24 5WE
Tel: +44(0)1954 233 100 Fax: +44(0)1954 233 101
Email: sales@camlab.co.uk Web: www.camlab.co.uk