



# COOLING FASCINATION ❄️

Refrigerators and freezers for professional use ❄️ ❄️

## Circulating cooler DLK 402

The smallest in its range. The DLK 402 circulating cooler serves as a universal cooler for the entire laboratory. A new method of casing construction provides for modern design, low weight and guarantees innovative technical features. Due in particular to the compact space-saving design with low-noise cooling unit, the freezer unit is suitable for direct use on the desk in the workplace.

Model	DLK 402
Temperature range	-10°C to +40°C
Control accuracy	+/- 1 K
<b>Cooling capacity in watts</b>	
at +30°C	380
at +10°C	270
<b>Circulation pump</b>	
Output	4l/min
Discharge head	1,5mWs
<b>External dimensions [WxDxH]</b>	26 x 37 x 40,5 cm
<b>Weight</b>	18,2 kg
<b>Ambient temperature</b>	+12°C bis +30°C
<b>Electrical connection</b>	230V / 50Hz
<b>Max. current</b>	1,5A
<b>Coolant tank</b>	4l



### Characteristics:

- Electronic temperature control with LED display
- Float switch for monitoring coolant level
- Level indicator on the front of the appliance
- Limit alarm when the permitted working area is exceeded
- Refrigeration unit: low noise, fully hermetically sealed, air-cooled, low maintenance
- Different hose connections included
- 3/4" filler neck and vent valve on the top of the appliance
- Casing made of rigid polyurethane foam
- Heat exchanger in non-ferrous metal version

### Accessories:

- Thermo-insulated hoses
  - Silicon Hoses: ID 10,0mm // to 0.5 bar // to -30°C // without fabric
  - Length on demand (maximum length: 10m)
- TP Thermal plate
  - For cooling and stabilising temperatures in conjunction with FRYKA DLK circulating coolers
  - two hose fittings made of stainless steel (suitable for DLK)
  - Stainless steel plates
  - Processing of refrigerated products during cooling possible
  - Low height: ergonomic work // can be placed on table // for installation
  - Separate installation of refrigeration unit and thermal plate possible
  - Plate size as desired
- HKF coolant 15.1
  - Recommended for all FRYKA circulating coolers // ready-made mixture // frost-proof up to -15°C

