




Speed, Performance, Value - The Right Balance

- Up to 12 Weighing units (including one custom unit)
- Date and time
- Full range tare
- Zero Tracking
- Multilingual display
- Large backlit display with dual text prompts
- Capacity tracker
- Large stainless steel pan
- Solid metal housing
- Overload protection
- Lock down mounting slot for Kensington™ type lock and cable
- Non-slip adjustable levelling feet
- Colour coded, sealed keypad
- Dual tare keys
- Removable draft shield supplied on 0.001g readabilities
- Modern low profile design
- RS-232 bi-directional interface
- GLP print outs
- External calibration
- Rechargeable battery
- AC Adapter
- Auto sleep / power down function to save battery life
- Weighing
- Parts counting
- Percentage weighing
- Check weighing
- Density determination
- Below balance weighing



Technical Specifications

Model	PGL 203	PGL 303	PGL 2002	PGL 3002	PGL 4001	PGL 6001	PGL 8001	PGL 10001	PGL 12001	PGL 15001	PGL 20001
Maximum capacity	200g	300g	2000g	3000g	4000g	6000g	8000g	10kg	12kg	15kg	20kg
Readability	0.001g	0.001g	0.01g	0.01g	0.1g	0.1g	0.1g	0.1g	0.1g	0.1g	0.1g
Repeatability (S.D.)	0.002g	0.002g	0.02g	0.02g	0.2g	0.2g	0.2g	0.2g	0.2g	0.2g	0.2g
Linearity (\pm)	0.004g	0.004g	0.04g	0.04g	0.4g	0.4g	0.4g	0.4g	0.4g	0.4g	0.4g
Pan size	145x125mm		192x192mm				400x300mm	192x192mm		400x300mm	
Net weight	4kg						6kg	4kg		6kg	

General Specifications

Typical stabilisation time	2-3 Sec
Units of measurement	Gram (g), Kilogram (kg), Carat (ct), Grain (GN), Newton (N), Pound (lb), Ounce (Oz), Pound-Ounce (lb-oz), Custom unit
Operating temperature	10°C 40°C
Power supply	AC adaptor and rechargeable battery pack
Calibration	External
Display	Backlit LCD display with dual digits and capacity tracker
Overall dimensions (wxdxh)	358x251x104mm

Key Accessories

- F1 5kg Calibration Weight for PGL 15001
- Density Determination Kit (mg models only)
- RS-232 to USB Interface Cable (must be ordered with 30740102)
- Security Lock and Cable
- Anti-Vibration Table
- Below Balance Hanger
- In-use Cover
- Dust Cover
- Adam Printer
- AdamDU - Data Collection Program
- RS-232 Cable
- Thermal Paper Roll for Adam Printer
- F1 1kg Calibration Weight for PGL 3002
- F1 1kg Calibration Weight for PGL 2002
- F1 200g Calibration Weight for PGL 303
- F1 2kg Calibration Weight for PGL 4001
- 2 x F1 5kg Calibration Weights for PGL 20001
- F1 200g Calibration Weight for PGL 203
- F1 5kg Calibration Weight for PGL 10001
- F1 5kg Calibration Weight for PGL 8001
- F1 5kg Calibration Weight for PGL 6001
- F1 5kg Calibration Weight for PGL 12001



www.camlab.co.uk
 24 Norman Way Industrial Estate, Over, Cambridge, CB24 5WE
 Tel: +44(0)1954 233 100 Fax: +44(0)1954 233 101
 Email: sales@camlab.co.uk Web: www.camlab.co.uk

www.WUa `UV'Ve'f _