

M SD 321	SD 321	SD 391	SD 481	SD 631	XL 391	XL 392	XL 481	XL 482	XL 483	XL 632	XL 633	XL 712	XL 713	XL 714
----------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------	--------

<b>Timer</b> Periodic reminder for the filter saturation test - 60 hours setting by default according to NFX 15 211: 2009 standard Range adjustable from 1 to 250 hours	N/A	S	S	S	S	S	S	S	S	S	S	S	S	S
<b>Flow monitor</b> Permanent ventilation control device according to EN 14175-2:2003 paragraph 8.2	N/A	S	S	S	S	S	S	S	S	S	S	S	S	S
<b>Air flow meter</b> Permanent air face velocity monitoring system	S	S	S	S	S	S	S	S	S	S	S	S	S	S
<b>Sampling port</b> for the filter saturation detection with color coded tubes (not supplied)	N/A	S	S	S	S	S	S	S	S	S	S	S	S	S
<b>Molecode S</b> Automatic alarm to detect filter saturation by solvents	N/A	O	O	O	O	O	O	O	O	O	O	O	O	O
<b>Work surface N°1</b> Metallic spill tray with tempered glass work surface	O	O	O	O	O	O	O	O	O	O	O	O	O	O
<b>Work surface N°2</b> In phenolic resin	O	O	O	O	O	O	O	O	O	O	O	O	O	O
<b>External lighting</b> IP 20 - 11 Watts - 400 Lux	O	O	O	O	O	O	O	O	O	O	O	O	O	O
<b>Internal lighting</b> IP 68 - 18 Watts - 500 Lux	N/A	O	O	O	O	O	O	O	O	O	O	O	O	O
<b>Benchcap</b> Fixed work bench	O	O	O	O	O	O	O	O	O	O	O	O	O	O
<b>Mobicap</b> Rolling cart	O	O	O	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
<b>Transparent back panel</b>	O	O	O	O	O	O	O	O	O	O	O	O	O	O
<b>Removable rear panel (maintenance access)</b>	N/A	N/A	N/A	N/A	N/A	O	O	O	O	O	O	O	O	O

S Standard      O Optional      N/A Non applicable