

New for 2009



Mobile Ductless fume hoods with modular filtration column

Expanded use - Modular design



Total operator safety

Compliant with the requirements of the ASHRAE 110:1995 & NF X 15-211:2009 standards

esp® program (Erlab Safety Program)*

100% energy savings

Very low energy consumption
Very easy to install, no ducting necessary

0 emission of toxic compounds

Toxic gases totally eliminated by filtration

* see back

One product - Many solutions



Distributed by

A single product to protect the user, the environment and your budget.

Coupon to return by mail or by fax to +44 (0)1954 233101

Please send me the following information :

captair® Flex™
Mobile Ductless fume hoods
with modular filtration column

captair® flow
Hepa filtered
enclosure.

captair pyramid®
Portable biological
investigation glove box.

captair® bio
PCR workstations.

captair® store
Ductless filtering
vented storage
cabinets.

asura® filters
Replacement filters
of all brands.

More information on www.camlab.co.uk/erlab

YOUR NAME _____
COMPANY NAME _____
ADDRESS _____
TEL _____
FAX _____
e-mail _____



985CC



New for 2009



Mobile Ductless fume hoods with modular filtration column

Expanded use Modular design

Flex™ Technology

A flexible/adaptable solution

A filtration technology allowing the use of liquids and solids individually or together with a single unit.

Enclosure dimensions from 80 cm to 1,80 m.

A cost cutting solution

Very low energy consumption.

Very easy to install and relocate, no ducting necessary.

Reduce carbon footprint, infrastructure, set up and running costs.

A safe solution

Compliant with the requirements of the following standards :

NF X 15 211:2009, ASHRAE 110:1995,...

An environment-friendly solution

No emission of toxic compounds into the atmosphere.



Erlab Safety Program

Essential for operator safety, our laboratory will assess the interactive behaviour of the molecules as well as the interactions with the filtration technology and will provide you with a **constant monitoring of the unit regarding the handlings performed.**



Call your ESP® specialist today and configure this unique solution to your requirements

Tel. : +44 (0)1954 233100
sales@camlab.co.uk



Distributed by

camlab

A single product to protect the user, the environment and your budget.



Worldwide specialist in lab filtration technologies since 1968. www.camlab.co.uk/erlab

FREEPOST
No Stamp
required

Camlab Limited

FREEPOST

Camlab House
Norman Way - Over
CAMBRIDGE
CB24 5YE

985CC

