



### Camlab Site Survey Form

Complete this page and any of the necessary attachments for the equipment being delivered and installed

Company:  Contact Name:

Delivery Address:  Contact Email:

Contact Number:

Equipment:

Quote Ref:

Avidity Science Rep

#### Delivery Access:

Door Dimensions  width  height

width  height

Lift Dimensions

Stairs How many steps?   
How many landings?  Size of landing   
Internal or External Staircase


Ceiling Height

Description/Drawing of Agreed Delivery Route:

Services Required




Socket for EACH water system & leak detector




Double socket required where system has seperate boosted storage tank (eg DUO systems)




Fused Spur



Isolator

PICO Only

Required PSU Cable Length 1000mm (POWE00002 - Standard)   
 2300mm (POWE00001 - extra cost)   
 Other (specify - extra cost)   
 Customer to supply extension



Water supply 3/4" Male BSP



Washing Machine Valve 3/4" male BSP

Other supply - please detail



Standpipe Drain



Sink / Tundish Drain



Floor Drain

Is the water system to be wall mounted? Yes  No

If yes to above, have you checked any plasterboard wall has been strengthened? Yes

Tanks

For completion if an Avidity Science tank (TANK30, TANK60, TANK100 or TANK100P-2) are to be installed

Wall mount Yes  No

Faucet Position Bottom  Middle  None

What will the system be supplying water to?   
 Specify type & quantity

Extra Requirements

Additional pre-treatment requirec Yes  No

Please provide details



**Additional Comments/ Drawings**