

# Freezer Boxes

## FREEZER STORAGE SOLUTIONS



200  
Place Boxes



50 Place  
Box

Height  
adjustable Boxes



100  
Place  
Boxes

## FOR ALL TUBES, MICROTUBES, VIALS AND CRYOVIALS

These freezer storage boxes are UK manufactured in 100% virgin materials. They offer perfect support for organised storage, transportation, freezing and autoclaving of your tubes.


Tower and tray racks are also available to facilitate easy handling.

For further details please contact us using the information overleaf.



### Features





- 50, 100, 200 place and height adjustable boxes in the range
- Manufactured in UK from 100% virgin polypropylene
- For use with vials, tubes and microtubes from 0.2 to 5ml
- Temperature range -80°C to 135°C
- Suitable for use in ultra-low freezers
- Suitable for autoclaving (15 minutes at 121°C)
- All are Stackable.
- All boxes also available with self adhesive grid which attaches to the underneath of the box
- Reinforced triple hinges prevent cracks or breaks
- Click-fit lock ensures box will not open accidentally
- Large clasp allows easy opening with gloves.



camlab

[www.camlab.co.uk](http://www.camlab.co.uk)

## ORDERING INFORMATION

	Description	Part Number	Colour	Pack Size
	50 Place Polypropylene Storage Box for 0.2 -1ml tubes up to 30mm in diameter.	RTP/72051-B	Blue	Pk of 5
		RTP/72051-G	Green	Pk of 5
		RTP/72051-N	Natural	Pk of 5
		RTP/72051-P	Pink	Pk of 5
		RTP/72051-Y	Yellow	Pk of 5
		RTP/72051-MX	Mixed Pack	Pk of 5
	100 Place Freezer Box for use with vials, tubes and microtubes from 0.5 to 2ml in volume and can hold tubes up to 53.5mm in diameter. Alphanumeric grid on lid for easy reference. Stackable, with hinged lid. All Black boxes available for use with light sensitive samples. A self adhesive grid is also available which attaches to the underneath of the box.  Dimensions: (HxWxD) 61x140x145mm	RTP/72101-B	Blue	Pk of 5
		RTP/72101-BK	Black	Pk of 5
		RTP/72101-G	Green	Pk of 5
		RTP/72101-N	Natural	Pk of 5
		RTP/72101-R	Red	Pk of 5
		RTP/72101-Y	Yellow	Pk of 5
		RTP/72101-MX	Mixed Pack	Pk of 5
		RTP/72101-FB	Neon Blue	Pk of 5
		RTP/72101-FG	Neon Green	Pk of 5
		RTP/72101-FO	Neon Orange	Pk of 5
		RTP/72101-FP	Neon Pink	Pk of 5
		RTP/72101-FY	Neon Yellow	Pk of 5
		RTP/72101-BI	Blue + Grid	Pk of 5
		RTP/72101-GI	Green + Grid	Pk of 5
		RTP/72101-NI	Natural + Grid	Pk of 5
		RTP/72101-RI	Red + Grid	Pk of 5
RTP/72101-YI	Yellow + Grid	Pk of 5		
RTP/72101-MXI	Mixed + Grid	Pk of 5		
	200 Place Polypropylene Storage Box for 0.5-2ml tubes. Alphanumeric grid on lid for easy reference. Stackable, with hinged lid. Natural box also available with self adhesive grid.	RTP/72200-N	Natural	Each
		RTP/72200-NI	Natural + Grid	Each
		RTP/72200-BK	Black	Each
	Height Adjustable Freezer Boxes with and without dividers. Height adjustable, uses 'stops' for heights of 55mm, 50mm, and 43mm. For use with vials, tubes and microtubes from 0.5 to 5ml volume. DIY grid available in same colours for you to choose how many divisions per box (from 4-81)	RTP/72700-B	Blue	Pk of 5
		RTP/72700-BK	Black	Pk of 5
		RTP/72700-G	Green	Pk of 5
		RTP/72700-O	Orange	Pk of 5
		RTP/72700-MX	Mixed	Pk of 5
		RTP/72700-N	Natural	Pk of 5
		RTP/72700-R	Red	Pk of 5
		RTP/72700-Y	Yellow	Pk of 5
		RTP/72701-DB	Blue dividers	Pk of 16
		RTP/72701-DBK	Black Dividers	Pk of 16
		RTP/72701-DG	Green Dividers	Pk of 16
		RTP/72701-DN	Natural Dividers	Pk of 16
		RTP/72701-DO	Orange Dividers	Pk of 16
		RTP/72701-DR	Red Dividers	Pk of 16
RTP/72701-DY	Yellow Dividers	Pk of 16		

## ALSO AVAILABLE FROM CAMLAB




Double Tamper Evident Bottle



Tubes® Tube Labels




Polycarbonate Cryo Boxes



Tamper Evident 1.5ml Tube



1.5ml Graduated Microtube



Camlab Limited, Camlab House, Norman Way, Over, Cambridge, CB24 5WE  
 Tel: +44 (0)1954 233110 • Fax: +44 (0)1954 233101  
 email: [plastics@camlab.co.uk](mailto:plastics@camlab.co.uk) • [www.camlab.co.uk](http://www.camlab.co.uk)