




## Polydisc GW and Polycap GW related products

Polydisc GW and Polycap GW have been developed for the preparation of larger volumes of groundwater samples for the analysis of dissolved heavy metals. The Polydisc GW and Polycap GW feature a stepped hose barb connection on the inlet and outlet for easy fitting with various tube sizes.



**Fig 3.** (A) Polydisc GW has a filtration area of 20.4 cm<sup>2</sup>. (B) Polycap GW is larger with a filtration area of 600 cm<sup>2</sup>.

## Ordering information

Product	Quantity	Code No.
<b>GD/XP Syringe Filters</b>		
GD/XP 25, 0.45 µm nylon	150	6970-2504
GD/XP 25, 0.45 µm nylon	1500	6971-2504
GD/XP 25, 0.45 µm PES	150	6994-2504
GD/XP 25, 0.45 µm PES	1500	6995-2504
GD/XP 25, 0.45 µm PVDF	150	6972-2504
GD/XP 25, 0.45 µm PVDF	1500	6973-2504
GD/XP 25, 0.45 µm PP	150	6978-2504
GD/XP 25, 0.45 µm PTFE	150	6974-2504
GD/XP 25, 0.45 µm DpPP	150	6993-2504

### Polydisc GW and Polycap GW

Polydisc GW, 0.45 µm nylon + quartz fiber	20	10 463 400
Polydisc GW, 0.45 µm nylon + quartz fiber	50	10 463 401
Polycap GW, 0.45 µm PES	1	6714-6004
Polycap GW, 0.45 µm PES	100	6724-6004

## cytiva.com/whatman

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. Whatman is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://cytiva.com/contact)

CY14462-01Oct20-DF



### Your Local Distributor



Camlab Ltd · 24 Norman Way Industrial Estate · Over · Cambridge · CB24 5WE · UK  
[sales@camlab.co.uk](mailto:sales@camlab.co.uk) · 01954 233 110 · [www.camlab.co.uk](http://www.camlab.co.uk)