



GW2045



GW2045 TECHNICAL FEATURES	All versions
Electronic control	Microprocessor
Standard programs stored	6
Customisable programs	2
Display with segments and warning light	yes
Clock and calendar	yes
Reprogrammable phases	8
Phase parameters	type of water, detergent quantity, target temperature, extension time in minutes, temperature
Tank internal temperature	from room temperature to 95°C
Accuracy	0,1°C
Temperature sensors in tank	PT 1000 CLASS B IEC 60751
Display	2 digits
Detergent dosing system	max. 1 powder dispenser and 2 liquid pumps
Detergent level sensor	optional
Alkaline powder detergent dispenser	yes
Pump to dispense acid neutralized liquid	yes
Pump to dispense alkaline liquid detergent	optional
Safety lock	yes, with electromagnetic release
Safety devices	safety thermostats, door interlock
Alarm display	30
Troubleshooting menu	yes, using PC connection
Programme editing	yes, using 2 spare programmes
AUXILIARY FUNCTIONS	
Traceability	storage of data for 100 most recently run programmes
RS232 serial port for PC connection	yes
RS232 serial port for printer	yes
Cycle storage	yes, WD-TRACE®/Excel format
Cycle file download	yes, WD-TRACE® format
WATER SUPPLY (PRESSURE 1,5-5 BAR)	
Cold water inlet	yes, max 42° F hardness
Demineralised water inlet	yes, < 20 µS/cm
Demineralised water booster pump	optional
Water softener incorporated	yes
Recirculation pump	200 l/min
WATER HEATING	
Electrical	3,0 kW max
STEAM CONDENSER	optional
DIMENSIONS LxDxH mm	
Outside (with built-in top)	450x620x850 (830)
Inside	380x480x590
Net weight (Kg)	61
STAINLESS STEEL	
Wash tank	AISI 316L
Exterior covering	AISI 304
ELECTRICAL POWER SUPPLY	
Max. voltage/power	1/N/PE 230V ~ 50Hz 3,3 kW
NOISE LEVEL	max 50 dB

