

GW6290



HIGH CAPACITY LABORATORY GLASSWARE WASHER (405 LITRES)
PROFESSIONAL LINE - LARGE DIMENSION GLASSWARE WASHERS
[Download PDF](#)

- ✓ washing up to 95°C - antibacterial thermal disinfection
- ✓ programs: 20 default + 20 custom
- ✓ multi-language TFT colour touch screen
- ✓ functioning protected by different levels password
- ✓ USB port for downloading washing cycles log
- ✓ peristaltic pump for liquid detergent
- ✓ peristaltic pump for liquid neutralizer
- ✓ possibility of installing up to 3 additional peristaltic pumps
- ✓ washing chamber and internal door made of AISI 316L stainless steel
- ✓ all internal piping made of AISI 304 stainless steel
- ✓ external panels made of AISI 304 stainless steel
- ✓ detergent cabinet made of AISI 304 stainless steel for 4 x 5 liters cans
- ✓ washing up to 3 levels
- ✓ external dimensions LxDxH: 900 x 800 x 1940 mm
- ✓ internal working dimensions LxDxH: 680 x 670 x 890 mm
- ✓ n. 2 independent washing pumps
- ✓ high efficiency drain pump
- ✓ built-in boiler for fast cycles preheating water
- ✓ cold water, hot water and demineralized water connections
- ✓ double electronic door lock system for user safety
- ✓ emergency opening system
- ✓ electrical connection: 3/N/PE 400 V - 50 Hz - 18 kW max
- ✓ 60 Hz version available
- ✓ without trolley

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350
Phone +39 0522 8211 - fax +39 0522 821453 - service.instruments@smeg.it

