

GW7015



HIGH CAPACITY DOUBLE SLIDING DOOR LABORATORY GLASSWARE WASHER (348 LT)
PROFESSIONAL LINE - DOUBLE DOOR GLASSWARE WASHERS
[Download PDF](#)

- ✓ washing up to 95°C - antibacterial thermal disinfection
- ✓ capacity up to 6 adjustable levels
- ✓ programs: 20 default + 20 custom
- ✓ multi-language TFT colour touch screen on both sides
- ✓ functioning protected by different levels password
- ✓ integrated built-in panel printer (optional)
- ✓ USB port for downloading washing cycles log
- ✓ n. 2 motors forced hot air drying system
- ✓ ethernet port for LAN connection (optional)
- ✓ n. 2 sliding down full glass door
- ✓ electronic door lock system for user safety
- ✓ high pressure hot air drying system
- ✓ high efficiency Ecoslim steam condenser (optional)
- ✓ peristaltic pump for liquid detergent
- ✓ peristaltic pump for liquid neutralizer
- ✓ possibility of installing up to 3 additional peristaltic pumps
- ✓ detergent cabinet with removable tray for storing up to 4 cans of 5 liters and one of 2 liters
- ✓ washing chamber made of AISI 316L stainless steel
- ✓ external panels made of AISI 304 stainless steel
- ✓ external dimensions LxDxH: 900 x 1000 x 1942 mm
- ✓ internal working dimensions LxDxH: 626 x 812 x 685 mm
- ✓ cold water, hot water and demineralized water connections
- ✓ n. 2 independent washing pumps
- ✓ floor drain connection
- ✓ drain pump (optional)

- ✓ electrical connection: 3/N/PE 380-400 V - 50 Hz - 20 kW max
- ✓ 60 Hz version available
- ✓ without racks

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350
Phone +39 0522 82111 - fax +39 0522 821453 - service.instruments@smeg.it

