

General Technical Data WTW Electrochemistry Instruments

	Digital			Analogue		
	Benchtop meters					
	inoLab®Multi IDS			inoLab®		
	Multi 9630	Multi 9620	Multi 9310	7320	7310	7110
Self test upon switching on	•	•	•	•	•	•
Splash Proof Housing	•	•	•	•	•	•
Display	Color graphic backlit	Color graphic backlit	Graphic backlit	Graphic backlit	Graphic backlit	LCD
Simultaneous Temperature Display	•	•	•	•	•	•
Temperature Compensation	•	•	•	•	•	•
TFK325 attachable (pH)	—	—	—	•	•	•
Membrane Keyboard with Defined PressurePoint	•	•	•	•	•	•
User Selectable Languages	•	•	•	•	•	—
Memory: DataSets	10000	10000	4500	5000	5000	—
Real Time Clock	•	•	•	•	•	—
GLPSupporting Functions	•	•	•	•	•	—
Identification No.	•	•	•	•	•	—
Sensor Serial Number, Manually	—	—	—	•	•	—
Digital Sensor Recognition	•	•	•	—	—	—
Calibration Record	•	•	•	•	•	—
Calibration Interval Selectable 1 ... 999 days	•	•	•	•	•	•
User Administration (via MultiLab®User)	•	•	•	—	—	—
USB-Interface/Mini-USB/A	• / • / • / • / — / —	• / • / • / • / — / —	• / — / • / — / • / —	• / — / • / — / • / —	• / — / • / — / • / —	• / — / • / —
PC Connection	•	•	•	•	•	—
Software MultiLab® Importer	•	•	•	•	•	—
Software Update via Internet	•	•	•	•	•	—
Firmware Update via Internet	•	•	•	•	•	—
Built-in Printer Option	—	—	•	•	•	—
Wireless ready	•	•	•	—	—	—
Sensor Validation	•, QSC	•, QSC	•, QSC	•	•	•
Dimensions without printer (H x W x D)	mm 240x280x70 (in. 9.45 x 11.02 x 3.15)	mm 240x280x70 (in. 9.45 x 11.02 x 3.15)	mm 190x240x80 (in. 7.48 x 9.45 x 3.15)	mm 190x240x80 (in. 7.48 x 9.45 x 3.15)	mm 190x240x80 (in. 7.48 x 9.45 x 3.15)	mm 190x240x80 (in. 7.48 x 9.45 x 3.15)
Dimensions with printer (H x W x D)	mm — (in. —)	mm — (in. —)	mm 190x290x80 (in. 7.48 x 11.42 x 3.15)	mm 190x290x80 (in. 7.48 x 11.42 x 3.15)	mm 190x290x80 (in. 7.48 x 11.42 x 3.15)	mm — (in. —)
Weight without printer	kg 2.5 (lb. 5.51)	kg 2.5 (lb. 5.51)	kg 0.8 (lb. 1.76)	kg 0.8 (lb. 1.76)	kg 0.8 (lb. 1.76)	kg 0.8 (lb. 1.76)
Weight with printer	kg — (lb. —)	kg — (lb. —)	kg 1.0 (lb. 2.20)	kg 1.0 (lb. 2.20)	kg 1.0 (lb. 2.20)	kg — (lb. —)
Universal Power Supply	•	•	•	•	•	•
Battery Operated (not with built-in printer)	—	—	•	•	•	•
Certificates	CE	CE	CE	CE	CE	CE
SET Available	•	•	•	•	•	•
Warranty	3 years	3 years	3 years	3 years	3 years	3 years

	Digital			Analogue			
	Portable meters						
	MultiLine IDS			Profiline			
	Multi 3630	Multi 3620	Multi 3510	Multi 3320	pH/Cond 3320	3310	3110
Self test upon switching on	•	•	•	•	•	•	•
Waterproof Housing	IP 67	IP 67	IP 67	IP 67	IP 67	IP 67	IP 67
Display	Color Graphic	Color Graphic	Color Graphic backlit	Color Graphic backlit	Color Graphic backlit	Color Graphic backlit	LCD
Simultaneous Temperature Display	•	•	•	•	•	•	•
Temperature Compensation	•	•	•	•	•	•	•
TFK325 attachable (pH)	—	—	—	—	—	•	•
Keypad/Acoustic prompts	•	•	•	•	•	•	•
User Selectable Languages	•	•	•	•	•	•	—
Memory: DataSets	10000	10000	4500	5000	5000	5000	—
Real Time Clock	•	•	•	•	•	•	—
GLPSupporting Functions	•	•	•	•	•	•	—
Identification No.	•	•	•	•	•	•	—
Calibration Record	10	10	10	•	•	10	—
Calibration Interval Selectable 1 ... 999 days	•	•	•	•	•	•	•
User Administration (via MultiLab®User)	•	•	•	—	—	—	—
USB Interface/Mini-USB/A	• / • / • / • / — / —	• / • / • / • / — / —	• / — / • / — / • / —	• / — / • / — / • / —	• / — / • / — / • / —	• / — / • / —	• / — / • / —
PC Connection	•	•	•	•	•	•	—
Software MultiLab® Importer	•	•	•	•	•	•	—
Firmware Update via Internet	•	•	•	•	•	•	—
Wireless ready	•	•	•	—	—	—	—
Sensor Validation	•, QSC	•, QSC	•, QSC	•	•	•	•
Dimensions (H x W x D)	mm 180x80x55 (in. 7.09 x 3.15 x 2.17)	mm 180x80x55 (in. 7.09 x 3.15 x 2.17)	mm 180x80x55 (in. 7.09 x 3.15 x 2.17)	mm 180x80x55 (in. 7.09 x 3.15 x 2.17)	mm 180x80x55 (in. 7.09 x 3.15 x 2.17)	mm 180x80x55 (in. 7.09 x 3.15 x 2.17)	mm 180x80x55 (in. 7.09 x 3.15 x 2.17)
Weight	kg 0.4 (lb. 0.88)	kg 0.4 (lb. 0.88)	kg 0.4 (lb. 0.88)	kg 0.4 (lb. 0.88)	kg 0.4 (lb. 0.88)	kg 0.4 (lb. 0.88)	kg 0.4 (lb. 0.88)
Universal Power Supply	•	•	USB, optional	USB, optional	USB, optional	USB, optional	—
Battery Operated/ Rechargeable Batteries	recharg. (NiMH)	recharg. (NiMH)	• / •	• / •	• / •	• / •	• / •
Certificates	CE	CE	CE	CE	CE	CE	CE
SET Available	•	•	•	•	•	•	•
Warranty	3 years	3 years	3 years	3 years	3 years	3 years	3 years



Camlab Limited

E: sales@camlab.co.uk

T: +44 1954 233110