

handvac handpipettor

Optimized controlling and flexible handling

Compact system for precise and comfortable operation

The ergonomic and flexible handheld pipettor can be used in combination with pasteur pipettes, plastic aspiration pipettes and a variety of other adapters needed for your application. A very precise control and pressuresensitive button on the top of the handvac controls the needed vacuum level. An accidental aspiration of cells is prevented.

Design Features

- Optimized Vacuum Level Controlling
- Handvac is fully autoclavable
- High quality components for trouble-free operation
- Suction regulation through press key on top of the HandVac
- Comfortable Vacuum Hose connection
- Wide accessory range




handvac system - complete tool set for aspiration

Suction regulation through press key on top of the HandVac

Comfortable Vacuum Hose connection

Perfect adaptation to human hand shape, Considers natural hand form with ergonomic features



Easy change of adapters and accessory


Working Areas



Aspirate



Filtrate



Rinse



Dispose


Accessories:



Rubber Adapter for Pasteur pipettes



1-Channel Adapter for disposable tips
Cat.No.: 828852-05



8-Channel Stainless Steel Adapter
Cat.No.: 828852-03



Table Top Stand for handvac and accessories



8-Channel Adapter for disposable tips
Cat.No.: 828852-04

WELCH-ILMVAC
APPLIED VACUUM TECHNOLOGY

Am Vogelherd 20, 98693 Ilmenau • Phone: +49 3677 604 0 • Email: welch-ilmvac@gardnerdenver.com

Package Name

Including



Aspiration System - biovac 106

Cat.No.: 112037

- biovac 106
- secretion glass receiver
2 liters, autoclavable
- handvac professional
vacuum hose



Aspiration System - biovac 106

Maxi Package

Cat.No.: 112037-02

- biovac 106
- secretion glass receiver
2 liters, autoclavable
- spare secretion glass
receiver, 2 liters
- handvac professional
on/off foot switch for
biovac
- vacuum hose



handvac professional

Cat. No.: 112580

- handvac handpipettor
complete
- 1-Channel Adapter for
disposable tips
- 8-Channel Adapter for
disposable tips
- 8-Channel Stainless Steel
Adapter
- rubber Adapter for
pasteur pipettes
- Table Top Stand
- vacuum Hose



handvac double

Cat.No.: 112580-02

- handvac handpipettor
- y - connector
- vacuum Hose



handvac basic

Cat.No.: 112580-01

- handvac handpipettor
- 1-Channel Adapter for
disposable tips
- rubber Adapter for pasteur
pipettes
- vacuum Hose