

D Series


NEW!

High Speed Micro Centrifuge

D3024



CE ISO 9001/13485
CERTIFIED


www.camlab.co.uk
camlab

Features

- ❑ User-friendly large LCD display provides all visual information
- ❑ Quiet and stable operation benefit lab environment
- ❑ Brushless motor drives quickly and effortlessly accelerates the rotor to set speed
- ❑ Fully autoclavable Bio-safe rotor
- ❑ Dual door interlock design, more safe
- ❑ Special double air cooling design to keep lower rotor temperature
- ❑ Unique temperature display to protect sensitive samples
- ❑ Gently braking at low speed, better separation

Specifications

Max. Speed	15000rpm (200-15000rpm) , increment: 10rpm
Max. RCF	21380×g, increment: 10×g
Rotor	2ml×24, 0.5ml×36, PCR8 serial tubes
Run Time	30sec~99min~HOLD (Continuous operation)
Driving Motor	Brushless motor
Safety Devices	Dual door interlock, Over-speed detector; Over-temperature detector; Automatic internal diagnosis
Power	Single-phase, AC 220V ~ 240V, 50Hz/60Hz, 5A AC 110V ~ 120V, 50Hz/60Hz, 5A
Noise level (dB)	58
Dimensions (mm)	(L) 280× (D) 364× (H) 266
Weight	12 kg
Additional Features	Speed/Acceleration switch; Short-time run function; sound-alert function

Ordering information

Cat No.	P/N	Description
912015130000	D3024	High Speed Micro-Centrifuge, with AS24-2 rotor kits, US plug, 110V/50Hz/60Hz
912115120000	D3024	High Speed Micro-Centrifuge, with AS24-2 rotor kits, Cn plug, 220V/50Hz/60Hz
912215120000	D3024	High Speed Micro-Centrifuge, with AS24-2 rotor kits, Euro plug, 220V/50Hz/60Hz
912315120000	D3024	High Speed Micro-Centrifuge, with AS24-2 rotor kits, UK plug, 220V/50Hz/60Hz


Accessories

19400002	AS24-2	Rotor kits with cover, fixing clips and O' ring, 15000rpm, 2ml*24, use with D3024 & D3024R
19400003	AS36-05	Rotor kits with cover, fixing clips and O' ring, 15000rpm, 0.5ml*36, use with D3024 & D3024R
19400004	AS4-PCR8	Rotor kits with cover, fixing clips and O' ring, 15000rpm, 4-PCR8, use with D3024 & D3024R
19500001	A02P2	0.2ml rotor adapter, use with A12-2 & AS24-2 rotors, 24pcs/pk
19500002	A05P2	0.5ml rotor adapter, use with A12-2 & AS24-2 rotors, 24pcs/pk
19500003	A02P05	0.2ml rotor adapter, use with AS36-05 rotor, 36pcs/pk

O.E.M. option also available

Available rotors and adapters

Rotor Cat. No.		Max. Speed (rpm)	Max. Rcf (×g)	Capacity	Bio-safe	Centrifuge	
						D3024	D3024R
AS24-2	Base info.	15000	21380	2ml×24	Yes	•	•
	Adapters	A02P2:0.2ml A05P2:0.5ml				•	•
AS36-05	Base info.	15000	21380	0.5ml×36	Yes	•	•
	Adapters	A02P05:0.2ml				•	•
AS4-PCR8	Base info.	15000	21380	4-PCR8	Yes	•	•



www.camlab.co.uk
camlab

Camlab Ltd
24 Norman way Industrial Estate
Over
Cambridgeshire
CB254 5wE

01954 233 110
www.camlab.co.uk
sales@camlab.co.uk