

KNF Lab Pump Accessories

	Pump Systems	KNF Code	Description
VP Series pump	VP1HV	032498	VP1HV Service Kit
	VP2	043241	VP2 Service Kit
	VP3	030687	VP3 Service Kit
	VP4	030771	VP4 Service Kit
	VP5	030771 x2	VP5 Service Kit. (2x 030771 required for VP5)
Mini Laboport and N86/N811 pumps	N86 KN.18	043241	Service Kit
	N86 KT.18	043242	Service Kit
	4ID Hose Barb	001936	Replacement 4mm ID Hose barb. Suitable for N86 pumps - Sold singly, 2 required per pump.
	N811 KN.18	044066	Service Kit
	N811 KT.18	044067	Service Kit
	PME0781	085544/000000	N86/N811 Footswitch
	6ID Hose barb	000360	Replacement 6mm ID Hose barb. Suitable for N811 pumps or as larger barbs for N86 pumps. - Sold singly, 2 required per pump.
	N816.3 KN.18	056315	Service Kit
	N816.3 KT.18	056316	Service Kit
	N816.1.2.KN.18	056315	Service Kit
	N816.1.2.KT.18	056316	Service Kit
	N838.3 KN.18	111952	Service Kit
N838.3 KT.18	111953	Service Kit	
Chemically Resistant Pumps	N810 FT.18	058077	Service Kit
	N820 FT.18	058078	Service Kit
	N840 FT.18	058079	Service Kit
	N810.3 FT.18	057357	Service Kit
	N820.3 FT.18	057358	Service Kit
	N840.3 FT.18	057359	Service Kit
	N842.3 FT.18	057359	Service Kit
Variable Speed Pumps	N920G interface cable	307757	N920G interface cable 2M, for use with VC900 controller
General purpose pumps and compressors	N022 AN.18	030305	Service Kit
	N022 AT.18	032490	Service Kit
	N022 AN.18	000350	Fine control valve & vacuum gauge
	N022 AN.18	000349	Fine control valve & pressure gauge
	N022 AN.18	000346	Filter/Silencer
	N026.3/1.2 AN.18	030313	Service Kit
	N026.3/1.2 AT.18	118953	Service Kit
	N035 AN.18	030307	Service Kit
	N035 AT.18	032489	Service Kit
	N035.3/1.2 AN.18	030315	Service Kit
	N035.3/1.2 AT.18	118974	Service Kit
	N035 range	000354	Fine control valve & vacuum gauge
	N035 range	000482	Fine control valve & pressure gauge
	N035 range	000352	Filter/Silencer
	N145 AN.18	030310	Service Kit
	N145 AT.18	032512	Service Kit
	N145.1.2 AN.18	030318	Service Kit
N145.1.2 AT.18	121293	Service Kit	
RC600/900 Accessories	SC920 G	306473/117578	Fully Automatic, speed regulated vacuum System with Wireless remote control & gas ballast - 20 l/min, <2 mbar, 100-240V, 50/60Hz - IP20
	SCC950	129656/117578	Fully Automatic, speed regulated Vacuum System with two wireless remote controls allowing independent vacuum supply for two rotary evaporators - 50 l/min, <2 mbar, 100-240V, 50/60Hz - IP20
	C900 Chiller	301230	C900 Chiller. Splashproof keypad with LED display. Temp range -10 to +40°C. Cooling capacity 250w max.
	VC900	304108/304181	Chemically resistant, portable digital vacuum controller. For use with any suitable vacuum source.
	VC900 Mounting Bracket	306221	Bracket for mounting VC900 on RC600/RC900 rotary evaporator tower.
	Protective Cover	127204	Protective cover for RC600/900 heated bath. With hinged access and auto alarm/stop function.
	Seal	113046	Seal RC600/900
	Coolant valve	300853	Coolant valve RC900 only
RC600/900 Glassware	Standard Condenser	128160	Standard condenser unit. Implosion protected
	Dry Ice Condenser	301696	Dry ice condenser unit for. Implosion protected
	Refill Valve vapour tube NS29/32	300639	Refill valve
	NS29/32	126059	Vapour tube NS29/32
	Collection Flasks	300557	100 ml. Implosion protected
		300558	250 ml. Implosion protected
		128158	500 ml. Implosion protected
		113939	1000 ml. Implosion protected
		113938	2000 ml. Implosion protected
	Evaporation Flasks NS29/32	113079	50 ml.
		113080	100 ml.
		113081	250 ml.
		113082	500 ml.
		128159	1000 ml.
		113083	2000 ml.
		113084	3000 ml.
	Powder Flask NS29/32	300588	500ml.
		300589	1000ml
		300590	2000ml
Foam Brake NS29/32	301114	Foam Brake	
Laboport Accessories	VC900*	304108/304181	New chemically resistant, digital portable vacuum control unit
	NP810/820	028129	Baseplate for N810/820
	NP840	028128	Baseplate for N840
	NR800	085482	Vapour Trap with tubing
	NR800	026225	Vapour Trap with no tubing
	NH800	026231	Condenser
	NG810	028476	Gas Balast for pump model N810
	NG820/840	028477	Gas Balast for pump model N820/N840
	NC800 A	046380/045115	Vacuum Controller with valves, cabling & electrical supply unit
	NC800 B	046381/045169	Vacuum Controller only with NC800 A with valves and cabling
	NC800A/B	046382	Double vacuum controller with valves, cabling and electrical supply unit
NZ800	051677/000000	Portable Vacuum Controller	
Simdos Accessories	FS60T Particulate filter	165210	Simdos 02 PVDF Liquid filter, Porocity 70um, Max flow rate 600ml/min, Connections UNF 1/4in - 28 Male/Female
	FS60X Particulate filter	165212	Simdos 02 PEEK Liquid filter, Porocity 35um, Max flow rate 600ml/min, Connections UNF 1/4in - 28 Male/Female
	FS25T Particulate filter	165211	Simdos 10 PVDF Liquid filter, Porocity 70um, Max flow rate 600ml/min, For tube ID 3.2/4mm.
	FS25X Particulate filter	165213	Simdos 10 PEEK Liquid filter, Porocity 35um, Max flow rate 600ml/min, For tube ID 3.2/4mm.
	Wall mounting plate	160473	Wall mount for Simdos 10 or 02
	Stand fitting	160474	Fitting to mount Simdos 1U or 02 on retort stand. (Retort stand not included)
	Footswitch	155872	Footswitch for remote start/stop of Simdos 1U or 02, RC or RCP versions only
	Pulsation damper	167817	For Simdos 10, KT versions
	Pulsation damper	167818	For Simdos 10, TT versions
	Pulsation damper	167819	For Simdos 10, FT versions
	Connection tube, Simdos 02	166335	Simdos 02 connection tube. 1/16" tubing with UNF 1/4-28, 1mtr for fine dosing