

## **PRODUCT INFORMATION**

**Model: DF3524L**

**Type: Labcold iShin Ultralow Temperature Freezer (-35 to -45°C)**

For laboratory and general use where storage at temperatures down to -45°C is required. This range offers low noise, a high level of construction, robustness and features coupled with great reliability and a comprehensive five year warranty just for good measure.



Shown with  
optional  
temperature chart  
recorder

**Capacity: 730 litres (25.5 cu.ft.)**

**Construction: Scratch resistant epoxy exterior; White enamel interior; (stainless steel optional)**

**Cabinet mounted on heavy duty castors**

**Inner compartments: 4; with insulated inner doors**

**Exterior dimensions (mm): 1980 high x 1125 wide x 1180 deep**

**Features:**

**Microprocessor temperature control with digital display**

**Audible and visible hi/lo temperature alarms with remote contacts**

**Lockable handle**

**Exceptionally low noise: approx 50dBa**

**Low frost build-up due to door design and frame heaters**

**Two access ports (1 rear, 1 side)**

**CFC and HCFC free**

**Temperature set-point range: -35 to -45°C**

**Maximum ambient temperature: 32°C**

**5 year full parts and labour warranty**

**Order reference: DF 3524L**