

## PRODUCT INFORMATION

**Model: DF9010**

**Type: Labcold ilShin Ultralow Temperature Freezer (-85°C)**

For laboratory and general use where storage at temperatures down to -85°C is required. This range offers low noise, a high level of construction, robustness and features coupled with great reliability and a comprehensive five year warranty just for good measure.



Shown with optional  
temperature chart  
recorder

**Capacity: 295 litres (10.3 cu.ft.)**

**Construction: Scratch resistant epoxy exterior; White enamel interior; (stainless steel optional)**

**Cabinet mounted on heavy duty castors**

**Inner compartments: 1; with insulated inner lids**

**Exterior dimensions (mm): 1031 high x 1590 wide x 790 deep**

### Features:

**Microprocessor temperature control with digital display**

**Audible and visible hi/lo temperature alarms with remote contacts**

**Lockable lid (with suitable padlock, not supplied)**

**Exceptionally low noise: approx 50dBa**

**Low frost build-up due to lid design and frame heaters**

**2 Access ports**

**CFC and HCFC free**

**Temperature set-point range: -50 to -85°C**

**Maximum ambient temperature: 32°C**

**5 year full parts and labour warranty**

**Order reference: DF 9010**