

PRODUCT INFORMATION

Model: RAFR05042 (White, mild steel)(GD suffix = Glass door)

RAFS05042 (Stainless steel)(GD suffix = Glass door)

Type: Labcold Advanced Laboratory Refrigerator

For laboratory and general use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.



Capacity: 129 litres (gross)

Construction: RAFR05042: White plastified mild steel interior and exterior.

RAFS05042: AISI 304 grade stainless steel interior and exterior.

GD suffix: Denotes refrigerator is fitted with glass door

Shelves supplied: 3; Plastic coated wire mesh 530 x 325mm

Exterior dimensions (mm): 850 high x 630 wide x 567 deep

Features:

Microprocessor temperature control with digital display to 1°C

Hi/Lo temperature alarm – audible

Fan circulated interior

CFC and HCFC free

Temperature set point range: +2 to +8°C

Automatic defrost with condensate removal

Maximum ambient temperature: +43°C (GD suffix = only +32°C)

Order reference: RAFR 05042, RAFS 05042

RAFR 05042GD, RAFS 05042GD